

H79969	Medeci System Plate - Ivory	H80005	Medeci System Plate - Blue
H80017	Medeci System Plate - Red	H69800	Medeci System Cup - Yellow
H70097	Medeci System Cup - Blue	H70103	Medeci System Cup - Red
H70115	Medeci System Cup - Ivory	F23938	Medeci Spout Top
F23914	Medeci Anti-spill Top		

⚠ General Safety Warning ⚠

Please make sure that you read and follow these instructions carefully. Keep these instructions in a safe place for future reference.

1. INTENDED USE

Medeci Systemware offers flexible solutions to a wide range of eating and drinking situations. Made from lightweight, high impact polymers, the bright contrasting colours make the range especially suitable for the visually impaired. The deep inner wall of the plate aids people who need to scoop up their food with one hand. The Medeci cup features two detachable, easy to hold, handles and a range of tops to reduce the risk of spillage and help people who have difficulty with conventional cups.



2. UNPACKING

Remove the packaging and check all parts carefully. Contact your supplier immediately if parts are missing or damaged and do not try and assemble.

3. IN USE

Medeci Plate

Medeci plates are designed to combine light weight and high strength and are appropriate for people who have difficulty eating from normal tableware.

The plate has a raised outer lip, one handed eaters will find this particularly useful as it can be used to maintain the food steady while loading it onto a fork or a spoon, when used in conjunction with a non-slip table mat (F17639, F17665 or L35324).


Medeci System Cup

Medeci cups are a lightweight, high strength cup designed to offer considerable choice. People who find drinking from normal tableware difficult will find the Medeci cup easier to use. The Medeci systemware cup range also offers the choice of two tops:

1. Anti-spill top (F23938)
2. Angled spout top (F23914).

3.1 Medeci System Cup

Each Medeci cup comes with two handles packed separately. They can be fitted into the grooves around the outside of the cup. The position of the two handles is a matter of choice for the user. One or both of the handles can be used; if both handles are used they can be positioned in different ways to conform to the abilities of the user's hands (see system cup handling options below).



3.2 Handle Removal

The handles can be removed from the groove by inverting the cup and tapping the handle sharply on a hard surface such as the edge of a table or work unit.

3.3 Medeci Cup Tops

1. Spout Top (F23938)

The spout top has an angled spout for people who find it difficult to cope with conventional cups or cups with anti spillage tops. The spout is angled at 45 degrees to the horizontal when the cup is partly tipped. At this angle the spout may be used to depress the tongue to assist in control during drinking.

2. Anti Spill Top (F23914)

The anti spill top is designed primarily for people with tremor or poor hand control. It has four holes drilled in the top, three in one side and one on the opposite side, the user may choose whether to drink through the side with three holes or the side with one hole according to the flow rate desired. The primary purpose of the anti spill top is to contain the drink within the cup. This is a benefit for people who suffer from a tremor or have limited hand control, and hence are at risk of spilling or splashing the contents of the cup.



4. CARE & MAINTENANCE

- Wash using normal domestic procedures
- Dishwasher safe*
- Care should be taken to ensure that surface temperatures of hot cabinets/plate warmers are below 1000C

- Although Medeci Systemware is extremely tough it is not scratch resistant and can, in normal use, be scratched relatively easily by the use of sharp cutlery.

*It is advised if cleaned in a dishwasher, they be placed on the top shelf and away from the heating element, and a lower temperature cycle selected.

Important:

- DO NOT Use in a conventional oven/grill
- DO NOT Use in a microwave oven/grill
- DO NOT Use abrasive cleaning or drying products.

5. PRODUCT DISPOSAL

This product must be thoroughly cleaned/sterilised prior to disposal.

This product is manufactured from plastic and can be disposed of through Local Authority Household Waste and Recycling Centres or through Local Trade Waste Collectors.

6. PRODUCT GUARANTEE

NRS Healthcare guarantee this product for a period of 12 months, from the date of purchase, against faulty materials and workmanship.

This guarantee does not cover; abuse, misuse, interference or tampering with any part of the product.

Any guarantee issued by NRS Healthcare does not affect your statutory rights.

NRS HEALTHCARE Sherwood House, Cartwright Way, Forest Business Park
Bardon Hill, Coalville, Leicestershire, LE67 1UB, UK.

Tel: Customer Services: 0345 121 8111

Email: customerservice@nrshealthcare.co.uk

NRS Healthcare is a leading supplier of community equipment and daily living aids.
Please visit one of our websites to view our full range of products:
www.healthcarepro.co.uk or www.completecareshop.co.uk