

USER INSTRUCTIONS

CORNER SHOWER STOOL

N05370

Non-Padded Seat

N05382

With Padded Seat

General Safety Warning

Please make sure that you read and follow these instructions carefully. Keep these instructions in a safe place for future reference.

Maximum User Weight



1. INTENDED USE

For use in a shower to provide a stable, seated support, where a user is unable to stand for long periods.

2. SAFETY

Warning: *The maximum user weight must not be more than 127 kg (20 st).*

Ensure an even floor surface to avoid rocking.

The three ferrules must be in contact with the floor at all times.

Care should be taken when leaning forward.

Make sure the seat is fully located on the frame.

Inspect at regular intervals for damage or excessive wear, if found or suspected, remove from use immediately and contact your supplier.

3. UNPACKING

Remove the packaging and check all parts carefully. Contact your supplier immediately if parts are missing or damaged and do not try and assemble.

4. HEIGHT ADJUSTMENT

Remove the plastic E-clips connecting the stool legs and leg extensions.

Slide the leg extensions to the required position.

Note: *Set all legs to the same height.*

Replace the plastic E-clips.

Note: *Make sure the stool is fully secure before use.*

5. PRODUCT CARE: CLEANING

Clean regularly, use a general purpose, neutral detergent and warm water.

Check the suitability of the cleaning product before use. Do not use abrasive cleaners or cloths as this will permanently damage the surface and invalidate any warranty.

Do not use any form of solvent as this will damage the surface and weaken the product and invalidate any warranty.

To disinfect, wipe with a dilute bleach solution, and rinse thoroughly with fresh water.

This shower stool is suitable for thermal decontamination and can be steam cleaned or autoclaved on cleaning cycles up to 80°C.

6. PRODUCT CARE: MAINTENANCE

Inspect the stool at regular intervals, ensure the product is stable and rigid.

Check all fixings are fully tightened and in good condition.

Check ferrules (feet) monthly for signs of wear.

7. PRODUCT GUARANTEE

NRS Healthcare guarantee this product for a period of 24 months, from the date of purchase, against faulty materials and workmanship.

This guarantee does not cover abuse, interference or tampering with any part of the products. Issues arising from fair wear & tear and customer misuse are not covered.

This guarantee covers domestic use only. If the product is used in commercial or semi-commercial settings such as care homes, schools, rental or communal accommodation etc, the guarantee is valid for one year only. Consumable parts are also guaranteed for one year only.

Any guarantee issued by NRS Healthcare does not affect your statutory rights.

NRS HEALTHCARE Sherwood House, Cartwright Way, Forest Business Park
Bardon Hill, Coalville, Leicestershire, LE67 1UB, UK.

Tel: Customer Services: +44 (0)345 121 8111

Email: customerservice@nrshealthcare.co.uk

NRS Healthcare is a leading supplier of community equipment and daily living aids.
Please visit one of our websites to view our full range of products:
www.healthcarepro.co.uk or www.completecareshop.co.uk