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
Samsoft Elite 180 patient lift

TABLE OF CONTENTS

WARNINGS	30
I - OVERVIEW	31
1.1 - Labelling & symbols used.....	31
1.2 - Intended use.....	32
1.3 - Service life.....	32
1.4 - Advance safety information.....	32
II - INSTALLATION	33
2.1 - General terminology.....	33
2.2 - Unpacking.....	33
2.3 - Assembly.....	34
III - FEATURES & ACTIVATION	35
3.1 - Overview of the electric system.....	35
3.1.1 - Screen and main interface.....	35
3.1.2 - Control box.....	35
3.1.3 - Remote control.....	36
3.2 - Starting up the SAMSOFT ELITE 180.....	36
3.3 - Raising/lowering the boom.....	36
3.4 - Opening/Closing the legs.....	37
3.4.1 - Manual version.....	37
3.4.2 - Electric version.....	37
3.5 - Performing an emergency stop.....	37
3.6 - Emergency manual lowering.....	37
3.7 - Warning signals.....	38
3.8 - Charging.....	39
3.9 - Battery.....	39
3.9.1 - General safety instructions for the battery.....	39
3.9.2 - Activating the battery.....	40
3.9.3 - Battery storage.....	40
3.10 - Activation.....	40
IV - PATIENT TRANSFER	40
4.1 - Before transfer.....	40
4.2 - Handling the patient lift.....	41
4.3 - Installing the slings.....	41
4.4 - Transferring the patient.....	42
V - COMPONENTS COMPATIBLE WITH THE SAMSOFT ELITE 180	42
5.1 - Range of compatible slings.....	42
5.2 - Available options.....	42
5.2.1 - Offset 4-point spreader-bar.....	42
5.2.2 - Wall charging station.....	42
VI - TRANSPORT & STORAGE	43
6.1 - Folding.....	43
6.2 - Vertical storage of the patient lift.....	44
6.3 - Separating the base from the rest of the patient lift.....	45
VII - MAINTENANCE	46
7.1 - Cleaning.....	46
7.2 - Disinfection.....	46
7.3 - Device inspection.....	46
7.4 - Sling inspection.....	46
VIII - SYSTEM SETTINGS	46
IX - STORAGE, SHIPPING & DISPOSAL	49
9.1 - Storage.....	49
9.2 - Shipping.....	50
9.3 - Disposal.....	50

X - WARRANTY CONDITIONS AND LIABILITY	50
XI - TECHNICAL SPECIFICATIONS	50
11.1 - General characteristics.....	50
11.2 - Component weights.....	51
11.3 - Dimensions.....	51
11.4 - Electric system.....	52
11.5 - Environmental conditions.....	52
11.6 - Operating forces for controls.....	52
11.7 - Electromagnetic compatibility.....	53
XII - BEFORE CALLING DRIVE DEVILBISS HEALTHCARE FRANCE CUSTOMER CARE SERVICE	55
XIII - MAINTENANCE BOOKLET	113

- **READ THESE OPERATING INSTRUCTIONS AND THE INSTRUCTIONS FOR ACCESSORIES CAREFULLY BEFORE USING THIS DEVICE.**
- **USE ONLY THE SLINGS AND ACCESSORIES SUPPLIED BY DRIVE DEVILBISS HEALTHCARE FRANCE.**
- **MAXIMUM PATIENT WEIGHT: 180 KG**
- **THIS DEVICE MUST ONLY BE USED BY QUALIFIED PERSONNEL TRAINED BY THE DISTRIBUTOR.**
- **NEVER USE A DEFECTIVE DEVICE.**
- **IT IS COMPULSORY TO INSPECT AND MAINTAIN THE DEVICE. THE INSPECTION AND MAINTENANCE MUST BE CARRIED OUT IN ACCORDANCE WITH §VII - MAINTENANCE, OR ACCORDING TO LOCAL RULES AND REGULATIONS.**

WARNINGS

It is the responsibility of the distributor to inform and train the user in the handling of the patient lift, detailing the precautions for use and the risks associated with inappropriate use.

Pinch points must be explained in detail. The distributor must ensure that the user fully understands the safety instructions. When assembling and using the device, we recommend that you keep your hands away from the pinch points.

In case of prolonged exposure of the device to a source of heat, sunlight through glass or a radiator, ensure that there is no risk of burns before using the device.

The device must be used in a sufficiently clear space so that the surroundings do not interfere with manoeuvring.

It is prohibited to use the device near flames or to smoke near the device.

The device must not be used outdoors. Do not use it in a damp environment and avoid splashing water. Never immerse the patient lift.

The assembly must be carried out by the distributor's specialists. The weight and dimensions of the components are such that there are risks during assembly and dismantling; it is therefore highly advisable to seek assistance and pay attention to the risk of pinching.

If you have any questions or require information about these operating instructions, the manufacturer's contact details can be found on page 30.

Any serious incident occurring in connection with this device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is based.

I - OVERVIEW

SAMSOFT ELITE 180 is available in two versions:

- (M) for the Manual version
- (E) for the Electric version




Fig. 1: Electric leg opening version on the left,
Manual leg opening version on the right

The **SAMSOFT ELITE 180** is characterised by its main function of lifting a patient, as well as by the opening of the rolling base and the possibility of being folded for transport or storage.

Unless stated otherwise, the technical specifications are the same.


1.1 - LABELLING AND SYMBOLS USED

SYMBOLS AND MEANING	
	Maximum user weight
	Consult the instructions for use
	Refer to the instructions for use
	Manufacturer
	Date of manufacture
	Catalogue reference
	Medical device
	Serial number
	UDI – Unique Device Identifier Unique manufacturing identification number
	European compliance - CE marking indicating compliance with MDR 2017/745
	Warning: the warnings in this manual highlight risks which, if ignored, could result in injury or death.
	Caution: The cautions in this manual highlight risks which, if ignored, could result in damage to or malfunction of the device.
	Warning: risk of crushing / pinching
	Recommendations and tips
	Temperature limit
	Humidity limit
	Atmospheric pressure limit
	Class II electrical protection
	Type BF applied part
	For indoor use only

	Waste electrical and electronic equipment (WEEE), at the end of its service life, this device must not be disposed of with household waste because it must be recycled in an approved centre.
	Triman – Product subject to recycling rules
	In use
	In storage

EXPLANATION OF THE SAFETY LABEL

1.2 - INTENDED USE

The mobile patient lift is a battery-powered device for lifting a person with reduced mobility, in a seated, semi-seated or semi-recumbent position, from one location to another (application area 12 36 03 according to ISO 9999 classification). Examples: seat, bed, toilet seat, lifting a patient from the floor. This device is mainly used in medical establishments (e.g. retirement homes, etc.) or at home and must only be used indoors on a flat, stable and smooth floor.

It must be operated by competent personnel trained in the use of this type of device (healthcare professionals or trained private individuals).

Body support is provided by slings, which are medical devices independent of the patient lift (see instructions for use of slings). Only those distributed by DDHF must be used.

The patient lift must be medically prescribed or approved in advance, in order to assess the patient's transferability and the user's ability to carry out a transfer.

1.3 - SERVICE LIFE

The patient lift has an expected service life of 10 years under normal conditions of use and subject to proper maintenance as recommended by Drive Devilbiss Healthcare France in this manual.

Actual service life may vary according to conditions of use and the replacement of certain parts.

Slings and devices compatible with the patient lift have their own service life, which is not described here (refer to their specific instructions).

1.4 - ADVANCE SAFETY INFORMATION

Do not use this product or any optional equipment without completely reading and understanding the instructions contained in this document and all other manuals associated with this device. Complete the training provided by the distributor on the use of this device correctly. If you have difficulty understanding the warnings and instructions, contact a healthcare professional, distributor or DDHF technical staff before attempting to use this equipment.

The information contained in this document may change without notice. Before use, make sure you have the most recent version of the document. The latest version of our documents is available on the Drive Devilbiss Healthcare website.

Check all the components of the device on receipt; if they are damaged, do not use the device. Contact your sales representative or Drive DeVilbiss Healthcare France directly for further instructions.

In the event of a ground pick-up, do not transfer the patient without medical advice. The patient's medical condition may have been worsened by the fall.

- The patient lift is for indoor use only
- Never charge the device in a damp environment or where water could splash onto it.
- If used in a bathroom, be sure to wipe off any water splashes and moisture with a soft cloth after use.
- Never immerse the SAMSOFT ELITE 180.

Do not use the SAMSOFT ELITE 180 beyond its maximum load of 180 kg.

- Pinch points of the body extremities are present in several places on the patient lift (Fig. 2).
- The spreader bar could move suddenly and cause injury.
 - Always keep body extremities well away from moving parts.
 - When the patient lift is positioned, pay attention to the position of the spreader bar in relation to the patient.


Fig. 2: Location of areas with pinch points on the SAMSOFT ELITE 180

There is a risk of pinching the feet between the wheel base and the floor. Pay attention to the position of your feet when handling the patient lift.

When handling the patient lift, pay constant attention to your own position, to that of the patient and any accompanying person, so as to avoid impact with a moving part.

Parts of the patient lift colliding with the floor or with obstacles (people, walls, furniture, etc.) can cause damage to the device and lead to injury.

- Never allow the device or any of its parts to hit the floor, a wall or any other obstacle.
- Always ensure that the remote control is properly stored when not in use.

Only accessories approved and distributed by Drive DeVilbiss Healthcare France: slings, scales, etc. may be used with the patient lift.

The **SAMSOFT ELITE 180** is only designed to transfer one patient at a time. Any use not described in these operating instructions is prohibited.

The patient lift and its compatible accessories may not have the same maximum permissible load. The value to be taken into account must always be the lowest of the maximum loads of the patient lift and its accessories. The value used must be reassessed as soon as an accessory is added or removed.

The **SAMSOFT ELITE 180** requires regular maintenance, the steps for which are described in this manual §VII - Maintenance.

If the emergency stop button malfunctions in a dangerous situation, remove the battery from the control box as a last resort to power down the system.

II - INSTALLATION

2.1 - GENERAL TERMINOLOGY

1	Mast	13	Safety latch trigger
2	Lifting beam	14	Magnetic mast locking shaft
3	Legs	15	Mast clamping handle
4	Spreader bar	16	Front wheels
5	Spreader bar hooks	17	Rear wheels
6	Boom	18	Emergency stop button
7	Lifting actuator	19	Grip handles
8	Remote control	20	Brakes
9	Control box	21	Opening actuator (electric version)
10	Battery		
11	Leg opening pedal		
12	Charger support		

Fig. 3: Parts list for **SAMSOFT ELITE 180**



Fig. 4: Part for **SAMSOFT ELITE 180** electric version:

2.2 - UNPACKING

Ensure that the various components of the appliance have not been damaged in transit. Ensure that the folded device and its accessories are present: battery, operating instructions.

- Unpack the device at floor level.
- Remove the box structure.
- Place the structure on level ground and lock the wheel brakes.
- Cut the clamp (tie) that keeps the boom folded on the mast.

Note: Our equipment is assembled and tested at our factory. It is then folded and packed in its box.


Fig. 5: Samsoft Elite 180 folded

2.3 - ASSEMBLY

Incorrect or incomplete assembly may result in injury or damage.

- Assembly should only be carried out by qualified personnel.
- Use only DDHF components for assembly.

- Remove the magnetic locking shaft from its hole in the top of the lifting beam.
- Straighten the mast until the automatic safety latch engages.


Fig. 6: 1) remove the magnetic locking shaft, 2) straighten the mast until the safety latch engages.

- Be careful to take account of the mast's own weight when rotating it.
- An automatic latch holds the mast in place until the magnetic locking shaft is replaced. Straighten the mast until you hear a metallic click indicating that the latch has locked. Hold the mast until the safety latch engages.

- Position the magnetic locking shaft in its hole at the base of the lifting beam. The magnetised part at the base of the magnetic locking shaft must come into contact with the surface of the lifting beam, the pin must be fully inserted...


Fig. 7: Inserting the magnetic locking shaft


Fig. 8: View of the correctly inserted magnetic locking shaft

- Fully charge the battery as described in §III.8 Charging.
- Switch on the device and raise the boom using the remote control, as described in §III.2 Switching on the patient lift. The spreader bar, which is connected to the mast by a magnet, will detach automatically. A releasing noise is heard during detachment and is completely natural.


Fig. 9: Spreader bar output direction

Make sure that the spreader bar and in particular its hooks are well clear of the cables near the mast during this phase of unfolding the boom to avoid any risk of the cables being caught and pulled out.

Do not unfold the boom near the patient. A sudden detachment of the spreader bar could hit the patient.

III - FEATURES & ACTIVATION

3.1 - OVERVIEW OF THE ELECTRIC SYSTEM


Fig. 10: Electrical assembly

3.1.1 - Screen and main interface

The **SAMSOFT ELITE 180** has a screen on the control box which allows you to monitor the actions in progress.


Fig. 11: Main interface

	Clock - displays hours and minutes. The clock built into the electronics keeps time for more than 3 days, even if the device's battery is removed.
	Lift cylinder cycle counter - counter based on the distance travelled by the cylinder, once the cylinder travels 300 mm (corresponds to the average distance travelled for a lift and lowering in a transfer) the counter increments by 1.

	Overload counter recorded on the lifting cylinder.
	Months in use counter - number of months in use / number of months defined between each maintenance reminder (defined in the settings).
	Battery level indicator - continuously variable. The level changes to orange when the battery is below 20%. The level changes to red when the battery is below 10%.
	Dynamic display of the current action of the lifting or leg opening functions (in case of electric leg opening) – with an animation showing which of the functions is currently in operation (in red on the animation).

3.1.2 - Control box

Emergency stop button:
The device has an emergency stop button above the control box. Never block or obstruct access to the emergency stop button!

In case of an emergency, the emergency stop button immediately cuts off all movements, functions and power to the patient lift.

The emergency stop button also serves to power down the system: the control box is no longer supplied with power, the screen goes blank and all movement stops. Press the emergency stop button to activate it when the lift is not in use, or during storage or transport.

Emergency stop button enabled

Emergency stop button disabled

To reset the emergency stop button, turn the button slightly clockwise until it is unlocked and pops back up.

Control box screen:

Boom lifting
(M1 cylinder)

Opening the legs
(M2 cylinder)

Boom lowering
(M1 cylinder)

Closing the legs
(M2 cylinder)

3.1.3 –Remote control

	Overload signal - orange light when the system is overloaded.
	Boom Lowering + raising
	Legs Closing + opening (Electric spacing version)
	Maintenance reminder signal - orange light when the time set in the settings between maintenance intervals is reached.
	Battery level indicator
	Energy saving mode (flashing green)
	Different battery level options - fully charged with 3 green bars lit and almost discharged with 1 green bar lit

All the buttons on the control box and remote control are controlled by a sustained action: If the button is pressed, the function is activated; if the button is released, the function stops. The buttons on the control box take precedence over those on the remote control.

The remote control has a magnet for easy storage. As soon as the remote control is not in use, it should be magnetically stored on the mast or handles. When doing so, take care to ensure that the cable does not get pinched between two parts or caught when moving.

3.2 - STARTING UP THE SAMSOFT ELITE 180

ACTIVATION PROCEDURE

- 1 Unlock the emergency stop button by turning it slightly clockwise
- 2 The system is powered up but is still off.
- 3 Press any button on the box or remote control to start the system.
- 4 The screen lights up and displays the Drive DeVilbiss logo for 3 seconds.
- 5 The main menu is displayed, and the device is ready to use.
- 6 After 2 minutes without interaction, the system enters power save mode: the screen goes blank, but the system is still on. The middle digit of the battery indicator on the remote control flashes to indicate energy-saving mode. As soon as a button is pressed, the energy-saving mode stops, and the electrical action starts immediately.
- 7 After 30 minutes without interaction, the system goes into full standby and switches off. To restart, press any button on the control box or on the remote control. The screen lights up, displays the Drive DeVilbiss logo for 3 seconds and then the main menu.

i Press the emergency stop button as soon as the patient lift is not in use to ensure that the system is completely powered down.

3.3 - RAISING/LOWERING THE BOOM

! If used incorrectly, the patient lift can tip over during lifting, endangering both the patient and staff. It is advisable to leave the rear wheel brakes disengaged to allow the device to assume its stable position at the start of the lifting.

To raise or lower the patient, use the respective buttons on the control box or remote control.

REMOTE CONTROL	CONTROL BOX	ACTION
		Press to lower the patient
		Press to lift the patient

i The patient lift moves electrically by pressing and holding the buttons; as soon as the buttons are released, the movement stops.

3.4 – OPENING/CLOSING THE LEGS

3.4.1 - Manual version

The manual opening version of the patient lift has a pedal-operated system.


Fig. 13: Overview of the leg opening pedal in closed and open position

- Press the left pedal **A** with your foot to open out the legs.
- Press the right pedal **B** with your foot to close the legs.

i Pressing the pedal all the way down opens/closes the legs completely. In these extreme positions, a system locks them in place.

3.4.2 - Electric version

The electric leg opening version of the patient lift has an independent electric cylinder for opening and closing the legs. To do this, use the respective buttons on the control box or remote control.

REMOTE CONTROL	CONTROL BOX	ACTION
		Press to open out the legs
		Press to close the legs

i The patient lift moves electrically by pressing and holding the buttons; as soon as the buttons are released, the movement stops.

- !** - Never rest your foot on the leg opening actuator to move the patient lift. Use only the handles.
- Make sure you use the electric leg opening on a floor free of obstacles.

3.5 - PERFORMING AN EMERGENCY STOP

In case of a dangerous situation, immediately press the emergency stop button on the top of the control box. The power supply will be interrupted, the screen will switch off, and all actions in progress will stop.


To restart the device, deactivate the emergency stop button by turning it clockwise. The power supply to the control box will be restored and the patient lift can be used normally again.

3.6 - EMERGENCY MANUAL LOWERING

! This system is an emergency device. Use it only in a dangerous situation, such as a suspended patient with the lift function inoperative.

In the event of partial or total loss of power, failure of the lift actuator or if the battery runs out in the middle of the transfer, the patient can be lowered manually. This is made possible by emergency manual lowering.


This device is located on the upper part of the lifting actuator rod and is indicated by a red ring.


lowering device

To lower the lifting actuator rod and lower the patient, turn the red ring clockwise.

3.7 - WARNING SIGNALS

Signal on the screen	Signal on the remote control	Meaning
 Surcharge		Overload signal: - As soon as an actuator (lifting or leg opening) is overloaded, the indicator light on the remote control lights up in orange and the display shows the overload signal. Release the actuator button, the indicator light goes out and the screen returns to the main interface. - If the overload occurs on the lifting actuator, the overload counter on the main interface will be incremented by +1. 
		Maintenance reminder signal: - When the maintenance reminder period is reached (set by the user in the settings, see §VIII- System settings), the indicator light on the remote control lights up in orange and the screen will display the maintenance reminder signal. The counters for the number of cycles and the number of months of use are displayed in orange. Maintenance should be carried out as soon as possible.   - Once the maintenance operation has been carried out and the maintenance reminder has been reset in the settings by the operator, the indicator light on the remote control will go out and the screen will return to the main interface.

Signal on the screen	Signal on the remote control	Meaning
	/	Reminder signal to unlock emergency stop while charging: Charging is only possible when the emergency stop button is deactivated. The signal disappears as soon as the user unlocks the emergency stop button.
	/	Battery not inserted signal: This signal appears if the control box is being charged with an uninserted or damaged battery. If, after checking, the battery is correctly inserted and the signal is still displayed, the battery is faulty and should be replaced immediately.
	/	Less than 20% battery signal When the battery reaches less than 20% capacity, this message appears on the screen for information purposes. The user must anticipate the need to charge the battery in the near future.
	/	Less than 10% battery signal When the battery reaches less than 10% capacity, this message will appear on the display and a sound signal will be heard while the actuators are being operated. Once this signal is displayed, it is still possible to raise and lower the lift. The user must charge the battery immediately once the transfer in progress has been completed.
	/	Overheating signal: - If the internal temperature of the control box exceeds 80°C, the display shows this signal and all movement stops. - Once the internal temperature of the control box drops below 45°C, the signal disappears and the patient lift can be used again.

	/	Opening actuator connection error signal: - When the lift actuator cable is mistakenly inserted into the socket provided for the spacing cylinder (M2), the display shows this error message to indicate that the connections should be checked. - It is not possible to operate the opening actuator while the error message is displayed.
	/	Lifting actuator connection error signal: When the lifting actuator button is pressed and its cable is incorrectly connected to the control box, the screen displays this error message. The user must check the connection of the lifting actuator cable (M1).

3.8 – CHARGING

There are two possible methods for charging the **SAMSOFT ELITE 180** battery: direct charging and indirect charging (with the wall-mounted charging station option, see §V.2.2).

i Note that it is impossible to operate the actuators under direct charging.


- Use only the charger supplied by Drive DeVilbiss Healthcare France.
- Never charge the device in a damp environment or where water could splash onto it.
- Unplug the charger when the device is fully charged and close the cover protecting the charging socket.

Direct charging procedure:


- Release the emergency stop button first.


- Open the rubber cover protecting the charging socket on the right side of the box and connect the charger connector.

- Connect the charger to a mains socket.
- The battery charging begins.
- The display on the control box shows the charging status, as illustrated below.

- When the charging is complete, the display will show that charging is complete. Unplug the charger from the mains socket. Disconnect the charger from the control box and close the protective rubber cover.

- The batteries must be charged every evening.
- Carry out this operation in a well-ventilated room.
- Check the condition of the charger and its cord before plugging it in.

3.9 - BATTERY

3.9.1 - General safety instructions for the battery

Stresses (dropping, damage, crushing, overheating, freezing) can alter the internal service life of the battery and cause short circuits in the long term, potentially long after the stress.

Batteries must therefore be stored and handled with care.

IMPORTANT:

If mishandled, the Lithium-Ion battery may explode, burn or emit toxic gases and cause chemical burns.

- Do not open, crush, dismantle or drop the batteries (no shocks).
- Keep out of reach of children.
- Do not expose batteries to heat or fire.
- Do not short-circuit a battery.
- Do not immerse batteries in water.
- In case of deep discharge, dispose of the battery (section IX – Storage, shipping, storage and disposal).

3.9.2 - Activating the battery

- Only use the charger supplied by Drive DeVilbiss Healthcare France.
- Charge the battery for at least 12 hours before using it for the first time.
- Make sure that the battery is securely locked when it is fitted in the control box; the two clips must be firmly engaged.

- The rechargeable battery only develops its full capacity after 5 to 10 charge/discharge cycles.
- Optimum battery performance is achieved at room temperature (20-25°C).
- The maximum operating temperature is +40°C.

3.9.3 - Battery storage

BATTERY STORAGE CONDITIONS	
• Storage temperature: -10°C to + 50°C	
• Store the device in a dry place away from dust and moisture. • Remove the battery and store it outside the control box. • The first charge must not take place more than 6 months after the manufacturing date shown on the label.	Charge the device at least every 6 months.

- Never leave lithium-ion batteries completely discharged.
- When the device is not in use for more than 2 weeks, remove the battery from the case.

3.10 - ACTIVATION

- Before using the device to transfer a patient, refer to the following sections of this manual:
 - §III Features & Activation
 - §I.4 Advance safety information
 - §IV Patient transfer

Before using the device, make sure that all the shafts and circlips are fully engaged in their respective holders.

The **SAMSOFT ELITE 180** is functional provided that a sling compatible with this device is used and that the maximum permissible load of the patient lift and its accessories is respected.

The **SAMSOFT ELITE 180** is designed to withstand splashes of water. The entire system is IPX6-rated, with the exception of the charger. The system can therefore be used in a bathroom, but the charger must be disconnected, removed from the patient lift and the rubber cover closed.

The **SAMSOFT ELITE 180** must be used exclusively indoors on a level, stable and anti-skid floor; it must be handled by personnel trained in the use of this type of device.

- Make sure that the patient lift is correctly assembled according to §II.3 - Assembly.
- Ensure that the actuator(s) respond correctly to commands from the remote control and/or the control box.
- Ensure that the batteries are sufficiently charged.
- The device is ready to use.

It is recommended to perform an off-load repetition cycle to get used to using the device, before transferring a patient for the first time.

- Press the emergency stop button as soon as the patient lift is not in use to ensure that the system is completely powered down.

IV – PATIENT TRANSFER

Actuators must only be used discontinuously: 2 minutes of operation and then 18 minutes of rest. Never exceed 2 minutes of continuous operation.

Any frequency of use exceeding the above cycle may result in damage to the lifting assembly (actuators & control box). This damage is excluded from the guarantee.

The magnetic locking shaft must always be inserted before operating the **SAMSOFT ELITE 180**.


Fig. 15: View of the correctly inserted magnetic locking shaft

4.1 – BEFORE TRANSFER

CHECKS BEFORE TRANSFERRING A PATIENT

1. Please read these instructions and the sling instructions carefully.
2. The patient's weight is compatible with the maximum load of this patient lift - 180 kg.
3. The device used to receive the patient (chair, bed, etc.) is compatible with the patient's weight.
4. The sling used is one of the slings compatible with this patient lift and is adapted to the patient (permissible load, restlessness, frailty, pathology, etc.).
5. The battery is charged and the device is fully operational.
6. The floor is flat, stable and unobstructed throughout the lift's journey.
7. The device brakes are functional.
8. Ensure that the sling is in good condition (see sling instructions).
9. The device is maintained in accordance with the instructions in this manual.
10. If the patient is restless, an attendant is on hand to help secure their movement.
11. The magnetic locking shaft is properly inserted and its magnet is in contact with the lifting beam.

4.2 – HANDLING THE PATIENT LIFT

Use of the legs:

- The legs can be placed in two positions: closed and open. In the case of the manual opening version, the mechanism includes automatic locking in the open or closed position.
- The closed position takes up minimal space and is recommended for movement. The opening of the legs is used to adapt to the environment and avoid obstacles (chassis, castors, foot-board, etc.). This makes it possible to get as close as possible to the patient during transfers.
- It is essential to pay attention to the position of the device in relation to the patient during transfer. The patient lift should be brought as close as possible to the patient so that the spreader bar is vertical to the patient's pelvis. If this is not done, there would be an imbalance when the device is lifted, which could lead to the device falling off and injuring the patient.

Using the handles:

- The handles must always be used to push/pull/orientate the patient lift.
- Never push/pull on the sling or any other part to move the patient lift.

Using the brakes:

- The brakes prevent any unwanted movement of the device. They should be used at the end of the transfer, when the patient is lowered.
- Never engage the brakes when lifting a patient. This would prevent the device from assuming its equilibrium position and create an imbalance during lifting which could cause the patient lift to fall.
- The patient lift brakes are also used for storage or during maintenance.

Improper use of this device may cause damage, injury or death.

The **SAMSOFT ELITE 180** is not a transport device. It is designed solely for transferring a patient from one surface to another (see §1.2 - Intended use).

- The patient lift must be used on a smooth, flat surface free of obstacles (door sill, etc.).
- Make sure you carry out the pre-transfer checks.
- Do not engage the lift's brakes when lifting.
- The spreader bar must be positioned vertically over the patient's pelvis. The legs can be adapted to the environment to get as close as possible to the patient.

4.3 – INSTALLING THE SLING

A sling that is not properly attached to the patient lift, or that is damaged or incorrectly adjusted, can cause the patient to fall or care staff to be injured.

- Read the instructions for your sling carefully and make sure you have been trained in its use.
- Only use a sling that is compatible with the **SAMSOFT ELITE 180** and approved by the prescribing doctor to ensure patient comfort and safety.
- Before each use, check the condition of the slings: there should be no loose seams, the fabric and slings should not be frayed, and there should be no tears or signs of wear (weft visible).

- Ensure that the patient's weight is compatible with the sling and the patient lift.
- Always ensure that the slings are securely fastened to the spreader bar before lifting the patient.
- When fitting the sling, make sure that it does not pinch the legs.

PREPARING THE SLING BEFORE TRANSFERRING A PATIENT

1. Mark the desired setting on the sling.
2. Attach the fasteners marked on the inside of the spreader bar hooks. The fasteners must be on the same setting on the left and right of the patient to ensure balance.
3. Make sure that the fasteners do not have any folds that could come undone suddenly during lifting and that they cannot come loose on their own from the hooks.
4. At the start of the transfer and before raising the patient completely, tension the sling by raising the patient until the slings are taut. Make sure that the patient is correctly positioned in the sling, that they are stable and that the fasteners are perfectly secure.

Installing the sling on a patient in bed:

- Position the patient on their side.
- Place the sling under the patient.
- Ensure that the dorsal part of the sling is correctly positioned (from the shoulders to the coccyx).
- Lay the patient on their back.
- Pass each leg support under each of the patient's thighs.
- Position the patient lift at the patient's level.
- Hook each sling buckle into the spreader bar hooks.
- The patient can now be transferred.

Installing the sling on a patient in a chair:

- Lean the patient forward holding them.
- Slide the sling down the patient's back up to the coccyx.
- Pass each leg support under each of the patient's thighs.
- Position the patient lift at the patient's level.
- Hook each sling loop into the spreader bar hooks.
- The patient can now be transferred.

Installation of the sling on a patient on the floor:

- Position the patient on their side.
- Place the sling under the patient.
- Ensure that the dorsal part of the sling is correctly positioned (from the shoulders to the coccyx).
- Lay the patient on their back.
- Pass each leg support under each of the patient's thighs.
- Position the patient lift at the patient's level so that the head and torso are inside the legs.
- Hook each sling loop into the spreader bar hooks.
- The patient can now be transferred.

In the event of a ground pick-up, do not transfer the patient without medical advice. The patient's medical condition may have been worsened by the fall.

4.4 - TRANSFERRING THE PATIENT

PATIENT TRANSFER

1. At the start of the transfer and before lifting the patient completely, tension the sling by lifting the patient until the slings are taut. Make sure that the patient is correctly positioned in the sling, that they are stable and that the fasteners are perfectly secure.
2. Lift the patient until their weight is fully supported by the lift and they can be moved away from their resting surface.
3. Using the handles, move the patient lift away from its initial position.
4. In the case of a transfer to a chair, bring the chair closer and position it under the suspended patient.
5. If the patient is being transferred to a bed, rotate the patient lift using the handles to position the patient above the bed.
6. Engage the lift's brakes and place the remote control in its holder.
7. Lower the patient to the new resting position. Make sure the patient is comfortable.
8. Engage the lift's brakes and place the remote control in its holder.
9. Remove the slings loops from the spreader bar hooks
10. Release the brakes.
11. Move the patient lift away.

In rare and exceptional cases, it is possible to use the patient lift to make a very short journey with the patient suspended (e.g. from bed to shower, from bed to commode chair, etc.). This movement must take place within the same room and the patient lift must not encounter any obstacles along the way (e.g. door sills, carpets, etc.). However, it is preferable to carry out this operation using a transfer chair.

To make these short movements with the patient suspended, use the handles to move and guide the patient lift, with the assistance of an attendant or helper who will secure the patient's movement. For the patient's comfort, we recommend lowering the boom as far as possible and turning the patient to face the lift operator. The patient's feet must rest on the lifting beam. As a result, the centre of gravity is at its lowest, which greatly improves the feeling of stability and safety.

- !** Moving a suspended patient is a high-risk task and requires great vigilance on the part of the user and/or carer. The patient could fall out of the sling or unbalance the device until it topples over.
- Do not push or pull on the slings or the spreader bar to move the patient lift.
 - Do not move an agitated patient.
 - Make sure to not collide with or roll over an obstacle (door sill, carpet, etc.).

V – COMPONENTS COMPATIBLE WITH THE SAMSOFT ELITE 180

The patient lift and its compatible accessories may not have the same maximum permissible load. The value to be taken into account must always be the lowest of the maximum loads of the patient lift and its accessories. It is necessary to reassess the value used as soon as an accessory is added or removed.

5.1 - RANGE OF COMPATIBLE SLINGS

Drive Devilbiss Healthcare France offers a range of slings compatible with the **SAMSOFT ELITE 180**.

Refer to our website for a presentation of this range:

Contact your local distributor to select the sling best suited to your needs.

5.2 - AVAILABLE OPTIONS

5.2.1 - Offset 4-point spreader bar

The standard spreader bar is 4-point in-line. A 4-point offset spreader bar is available as an accessory.


Fig. 16: Overview of the 4-point spreader bar

5.2.2 - Wall charging station

A wall-mounted charging station is available as an option for indirect charging: the battery is removed from the control box and connected to the charging station. This means that the patient lift can be used with another battery during the charging time.


Fig. 17: charging options: wall-mounted charging station on the left, table-top charging station on the right. A: Battery, B: Wall mounting, C: Charging station, D: Charger

Indirect charging procedure:

- Grasp the battery by its handle and press the locking clip.
- Separate the battery from the control box by pulling it towards you.

- Insert the battery into the charging station. The charging starts automatically.
- The charging station is fitted with a battery charge indicator light. Once the battery is charged, the indicator light changes from orange to green. The charging time is approximately 4 hours.

VI - TRANSPORT AND STORAGE

6.1 - FOLDING

The patient lift must be folded for transport and storage.

The principle of folding is similar to that of assembly (unfolding), but in reverse order:

- Press the down button on the remote control to fold the boom (see §IV.2 Electrical controls, for remote control operation). Hold the spreader bar with one hand to guide it until the spreader bar shaft makes contact with the position-holding magnet.


Fig. 18: Positioning the spreader bar on its magnet

- Lock the wheel brakes.
- Remove the magnetic locking shaft from its hole in the base of the lifting beam.


Fig. 19: Exit direction of the magnetic locking shaft

- Pull the safety latch trigger while holding the handles.


Fig. 20: Direction of actuation of the safety latch trigger, under the grab handles

- Using the grab handles and keeping the trigger pulled, rotate the mast until it touches the lifting beam.


Fig. 21: Pivoting mast and mast resting on the lifting beam

- Be careful to take account of the mast's own weight when rotating it. Hold the mast during rotation to prevent it from falling suddenly. Always adopt good ergonomic lifting practices.
- Risk of pinching.
- Be careful with cables to prevent them from being pinched between two parts.

- Position the magnetic locking shaft in its hole at the top of the lifting beam.


Fig. 22: Zoom in of the insertion of the magnetic locking shaft in its folded position hole

Before moving the device to the folded position, make sure that the magnetic locking shaft is properly inserted and that its magnet is in contact with the mast.

6.2 - VERTICAL STORAGE OF THE PATIENT LIFT

The patient lift can be transported and stored vertically. A procedure must be followed to ensure its stability.

The rear wheels are fitted with brakes that also lock in rotation. This function is used to create a stable surface between the two rear wheels and the bottom of the mast.

VERTICAL STORAGE PROCEDURE

1. Close the legs.
2. Ensure that the patient lift is in the folded position.
3. Rotate the wheels to position them in line with the legs, so that they are aligned in the cut-outs in the legs (see figure 24).
4. Lock the wheel brakes.
5. Lift the patient lift using the handles and tilt it until the base of the mast is in contact with the floor.


Fig. 23: Overview of the correct position of the rear wheel in the cut-out of the leg.


Fig. 24: Side view of patient lift tilted from folded to vertical position

- Be careful to take account of the mast's own weight when rotating it. Always adopt good ergonomic lifting practices.
- Beware of the risk of the patient lift tipping over.

6.3 - SEPARATING THE BASE FROM THE REST OF THE PATIENT LIFT

- Be careful to take account of the mast's own weight when dismantling it. Always adopt good ergonomic lifting practices.
- Beware of the risk of the patient lift tipping over.
- Always disconnect the cables from the control box and carry the charger and remote control separately when placing them in 2 parts. There is a risk of these components falling and/or the carrier tripping.

To make it easier to transport and handle the **SAMSOFT ELITE 180** at the point of use (e.g. climbing stairs), it can be dismantled into two parts.

- Put the patient lift in the folded position (§III.4 Folding). Place a protective cover (cloth, cardboard, etc.) under the boom.


Fig. 25: View of the folded patient lift and zoom in of the positioning of the magnetic locking shaft and the mast clamping handle

- Use a flathead screwdriver to open the cover at the base of the control box. This cover covers the cylinder cable and remote control connections.


Fig. 26: Opening the control box cover

- Disconnect the remote control, charger and cylinder cables from the control box. The remote control and charger should be carried separately in a pocket or bag to prevent them from falling out.
- Place a protective cover on the ground next to the base (cloth, cardboard, etc.) to receive the mast once it has been dismantled. The aim is to avoid marks on the boom or mast.
- Unscrew the mast clamping handle and remove it from its hole.


Fig. 26: Zoom in of the unfastened mast clamping handle

- Remove the magnetic locking shaft from its housing in the folded position.


Fig. 27: Zoom in of the magnetic locking shaft of the mast removed from its housing

- Separate the mast from the lifting beam and place the mast on the ground protection previously positioned.


Fig 28: View of the patient lift in 2 parts

- Re-position the magnetic locking shaft in a hole in the lifting beam.
- Re-position the clamping handle in its hole in the lifting beam.

Refer to the technical specifications for the mass of the two separate components.

Never lift the mast assembly by holding it by the cylinder or by the boom: the jack could be damaged and the magnet holding it in position could suddenly detach, causing the device or the user to fall. Lift the dismantled assembly by grasping the handles and the mast.

VII - MAINTENANCE

Unplug the mains cable and remove the battery before carrying out any work.

The **SAMSOFT ELITE 180** has been designed to minimise the need for daily maintenance, thanks to its robust and durable components. It is the user's responsibility to clean (see §VII.1 - Cleaning) and disinfect the device (see §VII.2 - Disinfection), and to carry out mandatory periodic inspections (see §VII.1 - Inspection of the device). These inspections must be recorded in the maintenance logbook.

Maintenance of the **SAMSOFT ELITE 180** should be carried out annually, or more frequently in the event of intensive use or in accordance with local regulatory requirements, by the technical teams. These must be staffed by technically trained personnel qualified in the maintenance of electromechanical equipment. A complete maintenance manual giving all the information required for preventive and corrective maintenance is available with your Drive DeVilbiss representative, on our website or via the following link:

7.1 – CLEANING

CLEANING DESCRIPTION	FREQUENCY
• Clean the appliance with a soft cloth soaked in a soap solution, followed by a thorough rinse and wipe. • Clean magnets of any metallic dust that may have attached to them. • Do not immerse the charger and charging station in water. • After cleaning, check the connection pins between the battery and the control box and wipe off any water that may have accumulated. • Do not use detergents that are too aggressive (e.g. pure bleach). • Do not flush with water.	Every month or whenever necessary

7.2 - DISINFECTION

DISINFECTION DESCRIPTION	FREQUENCY
• Disinfect the appliance using a detergent and disinfectant solution sprayed on its external surfaces.	After each use

The device is not designed to withstand automatic, high-pressure or steam cleaning and disinfection processes.

7.3 - DEVICE INSPECTION

INSPECTION DESCRIPTION	FREQUENCY
• Charge the battery every night	Every day
INSPECTION DESCRIPTION	FREQUENCY
• Place the device on level ground: ensure that all 4 wheels are supported. • Engage the brakes: it must not be possible to move the device (the floor must be flat and firm). • Run the machine off-load: ensure that the mast and legs move smoothly and quietly. • Use the mast handles: ensure that there is no excessive play between the mast and the lifting beam. • Check that the plastic covers and protective casings are not broken, cracked or missing.	Monthly or as soon as necessary
• Check the condition of the cables: mains cable, remote control, cylinder cables. If a cable shows the slightest sign of burning, wear or knocks, or if the wires are visible, immediately quarantine the device and have the repairs carried out by the technical team.	Monthly or as soon as necessary

If a control point is not validated during the inspection carried out by the user, the device must be taken out of use immediately. The user should contact their technical team, who will carry out a maintenance operation. Maintenance operations are detailed in the **SAMSOFT ELITE 180 maintenance manual.**

7.4 – SLING INSPECTION

A sling that is poorly attached to the patient lift, damaged or badly adjusted can cause the patient to fall or the user to be injured. Read the instructions for your sling carefully, make sure you have received training in the use of this device and follow the maintenance procedures in the sling instructions carefully.

DESCRIPTION OF SLING INSPECTION	FREQUENCY
• The slings must be inspected at least once a month, or more regularly in case of intensive use. • Ensure that the seams are not loose, that the web of the canvas and slings is not visible and that there are no torn parts.	Every month or whenever necessary



VIII - SYSTEM SETTINGS

A number of menus can be accessed to adjust the settings of the patient lift. These menus can be accessed using a combination.

The menus can only be accessed using the buttons on the control box.

NAVIGATING THE MENUS	
 	Direction arrows: move forward / backward one menu at a time
	OK: Enter highlighted menu / Confirm action
	RETURN: Exit selected menu / Cancel
COMBINATION OF ACCESSES TO THE MENU LIST	
Simultaneously press  and  for 3 seconds	

SETTINGS, USER & MAINTENANCE

LANGUAGE SETTINGS	 Menu overview from the menu list  Current menu overview	Different languages can be selected for the device: French, English, German and Chinese. • The currently selected language is shown framed by two grey lines in the menu list overview. • Press OK to enter the highlighted menu. The list of available languages appears blurred. • Use the arrow keys to select the desired language. • Press OK to select a language: the interface text changes and you exit the menu. Press BACK to retain the current language and exit the menu.
	 Menu overview from the menu list  Current menu overview	The patient lift clock can be set from year to minute. • The date currently used by the device is shown from the preview in the menu list. • Press OK to enter the highlighted menu. The pre-selected setting is highlighted in blue. The default setting is year. • Select the setting to be changed using the arrow keys. • Press OK to select and modify it. Use the arrow keys to select the desired value. If you hold down a direction arrow, the values increase/decrease rapidly. • Press OK to save the new setting value and return to the menu. Press RETURN to cancel the change and return to the menu. • Press BACK to return to the menu list.

Settings, user & maintenance																																		
HISTORY OF USE	USAGE HISTORY <table><tr><th></th><th>Month</th><th>Cycle</th><th>Overload</th></tr><tr><td>1</td><td>-</td><td>-</td><td>-</td></tr><tr><td>2</td><td>-</td><td>-</td><td>-</td></tr><tr><td>3</td><td>-</td><td>-</td><td>-</td></tr></table> Menu overview from the menu list USAGE HISTORY <table><tr><th></th><th>Month</th><th>Cycle</th><th>Overload</th></tr><tr><td>4</td><td>-</td><td>-</td><td>-</td></tr><tr><td>5</td><td>-</td><td>-</td><td>-</td></tr><tr><td>6</td><td>-</td><td>-</td><td>-</td></tr></table> Current menu overview		Month	Cycle	Overload	1	-	-	-	2	-	-	-	3	-	-	-		Month	Cycle	Overload	4	-	-	-	5	-	-	-	6	-	-	-	This menu shows the duration of use in months, the number of cycles completed and the number of overloads recorded. A new line is created for each maintenance reminder reset (see MAINTENANCE REMINDER RESET menu). Up to 9 lines can be stored in memory, beyond which the 1st line will be overwritten. • The first 3 lines are shown from the menu list. • Press OK to enter the highlighted menu. • Navigate between the different lines using the arrow keys. • Press BACK to return to the menu list.
		Month	Cycle	Overload																														
1	-	-	-																															
2	-	-	-																															
3	-	-	-																															
	Month	Cycle	Overload																															
4	-	-	-																															
5	-	-	-																															
6	-	-	-																															
MAINTENANCE REMINDER RESET	MAINTENANCE REMINDER RESET After performing maintenance Menu overview from the menu list MAINTENANCE REMINDER RESET After performing maintenance Press OK to confirm Reset overview before confirmation MAINTENANCE REMINDER RESET Overview of completed reset	A maintenance reminder can be activated (see MAINTENANCE REMINDER. CONFIG menu). Once the set limit has been reached, a maintenance reminder signal is displayed on the screen and remote control, informing the user of the need to carry out scheduled maintenance. Once this has been done, this menu can be used to start a new cycle and record a new line in the use history. • Press OK to enter the highlighted menu. • Press OK as shown to initiate the reset. • When is displayed, the reset is complete. The main interface counters are reset to zero. • Press BACK to return to the menu list.																																

Settings, user & maintenance

MAINTENANCE REMINDER CONFIGURATION

MAINTENANCE REMINDER SET.

Month

^

-

v

Cycle

^

-

v

Menu overview from the menu list

MAINTENANCE REMINDER SET.

Month

^

-

v

Cycle

^

-

v

Current menu overview

The maintenance reminder can be triggered by two cursors: duration in months or number of cycles carried out. If both cursors are used, the first one reached triggers the maintenance reminder display. The duration in months between two maintenance sessions is defined by the month setting. The number of cycles between two maintenance operations is defined by the cycle setting. The value of a cycle is explained in §III.1.1 - Screen and main interface.

- Press **OK** to enter the highlighted menu. The pre-selected setting is highlighted in **blue**. The default setting is month.
- Use the **arrow keys** to select the setting to be modified.
- Press **OK** to select and modify it. Use the arrow keys to select the desired value. If you hold down a direction arrow, the values increase/decrease rapidly.
- Press **OK** to save the new setting value. Press **RETURN** to cancel the change.
- Press **BACK** to return to the menu list.
- By default, the month and cycle settings are set to "-". They are therefore inactive until the user selects another value. It is possible to select a value for one setting and leave the other inactive with a value of "-".

The tools presented above, Maintenance History & Maintenance Reminder, should in no way replace standard maintenance monitoring documents and procedures.

- The Maintenance History does not replace the maintenance tracking document.
- The Maintenance Reminder does not provide a self-diagnosis of the device's need for maintenance.
- It is the user's responsibility to follow the maintenance rules in these instructions or in local regulations. These functions are simply tools to help the user in this task.

IX - STORAGE, SHIPPING AND DISPOSAL:

9.1 - STORAGE

STORAGE CONDITIONS	FREQUENCY
Before storing the device: • Clean the appliance thoroughly. • Ensure that the device is working properly.	Before each storage
Storage conditions: • Store the device in a dry place away from dust and moisture. • Storage temperature: -10°C to + 50°C • Humidity: 20% to 80%. • Atmospheric pressure: 860 hPA to 1060 hPA • Altitude: up to 2000m	For each storage

STORAGE CONDITIONS	FREQUENCY
Battery storage conditions: • Remove the battery to store it outside the control box. • The first charge must not take place more than 6 months after the manufacturing date shown on the label.	Charge the device at least every 6 months
Before activation after storage: • Charge the battery. • Ensure that the device is working properly	For each storage

9.2 - SHIPPING

TRANSPORT CONDITIONS

- The appliance must be returned or transported in its original packaging.
- The lift's battery is UN38.3 certified. The standard professional practices for transporting the device must be put in place.

9.3 - DISPOSAL

DISPOSAL CONDITIONS

- Defective parts, packaging and battery packs must be disposed of at approved waste disposal centres.
- The device must not be disposed of with household waste, but returned to an appropriate collection point or to distributor. The device is subject to Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and placed on the market after 13 August 2015.
- By following these instructions, you are doing your bit for the environment, helping to preserve natural resources and protect human health.

X - WARRANTY CONDITIONS AND LIABILITY

The **SAMSOFT ELITE 180** is under warranty for 5 years for the structure and cylinders, under normal use in accordance with the instructions in this manual. The warranty does not cover negligence or handling errors. The remote control and control box are under warranty for 2 years, and the battery pack for 1 year.

Any modification or addition of accessories without the written authorisation of Drive DeVilbiss Healthcare France will de facto exclude the device from the warranty and from Drive DeVilbiss Healthcare France's liability. Furthermore, we cannot be held liable for the use of slings not distributed by Drive DeVilbiss Healthcare France.

- The mandatory periodic inspections (§VII.3 - Inspection of the device) must be carried out and recorded in the patient lift maintenance logbook.
- Annual maintenance, or more frequent maintenance in the case of intensive use or in accordance with local regulations (§Maintenance manual), must be carried out by the technical teams and recorded in the patient lift's maintenance logbook.

Outside the warranty period, we cannot be held liable for any failure to carry out these inspections, servicing and maintenance. Only spare parts and components supplied by Drive DeVilbiss Healthcare France may be used for repairs.

- Any serious incident occurring in connection with this device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is based.

XI – TECHNICAL SPECIFICATIONS


11.1 - GENERAL CHARACTERISTICS

Maximum working load (patient + sling)	180 kg
Electrical operation	On battery
Type of battery charge	Standard: Internal, directly on device Option: External after removal of the battery pack, and connection to a charging station sold as an accessory
Battery charge duration	< 4h
Capacity of the device with 85 kg load, new and charged battery	More than 40 full cycles of lifting/lowering
Off-load lifting time	33 s
Off-load lifting speed	39 mm/s
Lifting time at max. working load	42 s
Lifting speed at max. working load	31 mm/s
Time for electrical opening of legs at maximum working load	(E): 4 s
Service life in number of lifting/lowering cycles at maximum working load	10000 cycles
Structure materials	Steel
Structure cladding	Epoxy powder paint


11.2 - COMPONENT WEIGHTS

Weight of wheel base	(M): 18.9 kg (E): 19.2 kg
Weight of mast + boom without spreader bar	20.9 kg
Total weight without spreader bar	(M): 39.8 kg (E): 40.1 kg
Weight of 4-point standard spreader bar in line	1.7 kg
Mast weight + boom with 4-point standard spreader bar in line	22.6 kg
Total weight with 4-point standard spreader bar in line	(M): 41.5 kg (E): 41.8 kg
Weight of optional offset 4-point spreader bar	2.6 kg

11.3 - DIMENSIONS

Unfolded OVERALL dimensions (in use) (W x D x H)	1325 ^{+/-10} x 650 ^{+30/0} x 1365 ^{+/-10} mm
Folder OVERALL dimensions (transport or storage) (W x D x H)	1260 ^{+/-10} x 650 ^{+30/0} x 425 ^{+10/-35} mm
EXTERNAL dimensions of packaging box (LxPxH)	1295 x 675 x 450 mm


A	Maximum range at 600 mm (reference height)	594 mm ±20
B	Maximum range from the base	608 mm ±20
C	Range from base with legs 700 mm apart	(M): 270 mm ±10 (E): 208 mm ±10
D	Minimum distance (open legs) from the wall to the CSP (central suspension point) at the maximum range	(M): 270 mm ±20 (E): 251 mm ±20
H	Minimum distance (open legs) from the wall to the CSP at the maximum height	(M): 532 mm ±20 (E): 513 mm ±20
T	Minimum distance (open legs) from the wall to the CSP at the minimum height	(M): 555 mm ±20 (E): 535 mm ±20
K	Maximum height of the CSP	1982 mm ±20
L	Minimum height of the CSP	672 mm
M	Lifting stroke	1310 mm ±20
K1	Maximum height of connection points	1782 mm ±20
L1	Minimum height of connection points	472 mm
N	Base height	119 mm ±1
Y	Ground clearance	(M): 27 mm ±1 (E): 26 mm ±1
P	Maximum internal width	(M): 1020 mm ±10/-40 (E): 1120 mm ±30/-40
Q	Internal width at maximum range	(M): 910 mm ±10/-30 (E): 991 mm ±25/-30
R	Minimum internal width	488 mm ±30/0
S	Front wheel diameter (80 mm wide, without brake)	Ø100 mm

T	Rear wheel diameter (80 mm wide, with brake)	Ø 100 mm
Y	Turning diameter (closed legs)	Ø 1430 mm
Z	Maximum leg opening angle	(M): 33° +1°/-2° (E): 40° +2°/-3°

General dimensional tolerances: +/- 5%.

NOTE: Dimensions shown with device unloaded

11.4 - ELECTRIC SYSTEM

Power supply	By dedicated battery pack
Battery type	Lithium-ion
Battery voltage	25.2V DC
Battery capacity	3.25 Ah / 81.9 Wh
Maximum load voltage, frequency and current	100 - 240V AC, 50 - 60Hz, Max 0.8A
Box output current	Max 13A
Duty factor	10%, 2min ON / 18min OFF @25°C
Protections	- Against battery overheating - Against over currents / overloads
Safety features	- Emergency stop button - Emergency manual lowering (in the event of a remote control fault) - Factory-set current limiter (to protect against overloads)
Charger input for internal battery charging	Yes

Protection class	Cylinders, box, remote control battery (in use): IPX6 Charger (charging): IP20
Insulation class	Class II equipment
	Type BF applied part Applied part compliant with the requirements for protection against electric shock of IEC60601-1
Noise level	≤55 dB(A)

11.5 - ENVIRONMENTAL CONDITIONS

	OPERATING CONDITIONS	STORAGE CONDITIONS
Temperature	From + 5°C to + 40°C	From -10°C to + 50°C
Humidity	From 20% to 80%, no condensation	
Atmospheric pressure	From 860 hPa to 1060 hPa	
Altitude	Up to 2000 m	

11.6 - OPERATING FORCES FOR CONTROLS

TYPE OF ACTION	FUNCTION OF THE CONTROL	EFFORT REQUIRED
With one finger	Control box and remote control buttons	< 5 N
With one hand	Emergency stop button	< 105 N
Rotary control	Emergency stop button (for emergency stop release)	< 1.9 Nm
	Ring on lifting cylinder for emergency manual lowering	< 1.9 Nm
With one foot	(M): Foot pedal for manual leg opening	< 300 N
	Brake pedal on castors	< 300 N

11.7 - ELECTROMAGNETIC COMPATIBILITY

In order to ensure the electromagnetic compatibility requirements to prevent all non-hazardous situations of our product, the EMC standard IEC 60601-1-2 - AMD1 has been verified.

The use of this equipment attached to or stacked with other equipment must be avoided as this could result in operating errors. If such use is absolutely necessary, this equipment and those around it must be monitored to ensure that they are operating normally.

The use of accessories, transducers or cables other than those supplied may cause an increase in electromagnetic emissions or a reduction in the electromagnetic immunity of the **SAMSOFT ELITE 180**, and cause problems during use.

Portable RF (radio frequency) communication equipment must not be used within 30 cm (12 inches) of any part of the power supply, cables or charger. Otherwise, the performance could be adversely affected.

Manufacturer's declaration and directives:

The **SAMSOFT ELITE 180** is suitable for use in the following environment. The user must ensure that it is only used in an appropriate electromagnetic environment, as specified.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - DIRECTIVES
RF emissions - CISPR 11	Group 1	The SAMSOFT ELITE 180 uses RF energy solely for its internal operation. As a result, its RF emissions are very low and are unlikely to cause interference with the surrounding electronic environment.
RF emissions - CISPR 11	Class B	The SAMSOFT ELITE 180 is suitable for use in all establishments, including domestic ones and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.
Harmonic current emissions - IEC 61000-3-2	Class A	
Voltage change, voltage fluctuation and flicker emissions - IEC 61000-3-3		

IMMUNITY TEST	IEC 60601-1-2 ed.4.1 COMPLIANCE	LEVEL OF COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - DIRECTIVES
Electrostatic discharge (ESD), IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be made of wood, concrete or ceramic tiles. If the floors are covered with a synthetic material, the relative humidity must be at least 30%.
Radiated RF electromagnetic field, IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	Portable and mobile RF communications equipment must not be used in close proximity to any part of the SAMSOFT ELITE 180 , including cables.
Proximity field immunity to RF wireless communication equipment	Level in the table below	Level in the table below	
RF line disruption, IEC 61000-4-6	3 V 0.15-80 MHz 80% AM at 1 kHz 6 V in ISM and amateur radio bands [0.15 MHz - 80 MHz].	3 V 0.15-80 MHz 80% AM at 1 kHz 6 V in ISM and amateur radio bands [0.15 MHz - 80 MHz].	

IMMUNITY TEST	IEC 60601-1-2 ed.4.1 COMPLIANCE	LEVEL OF COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - DIRECTIVES
Rapid transient / electro-static burst, IEC 61000-4-4	±2 kV for AC supply lines ±1 kV for input and output lines	±2 kV for AC power supply lines ±1 kV for input and output lines	The quality of the power supply must be the same as in a home, business or hospital environment.
Overvoltage, IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line ±2 kV line to earth	
Magnetic field at mains frequency, IEC 61000-4-8	30 A/m (50-60 Hz)	30 A/m (50-60 Hz)	The magnetic field at the frequency of the network must be the same as that of a home, business or hospital environment.
Voltage dips, short interruptions and voltage variations on the power supply, IEC 61000-4-11	0% UT for 0.5 cycle in increments of: 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT for 1 cycle and 70% UT for 25/30 cycles at 0°.	0% UT for 0.5 cycle in increments of: 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT for 1 cycle and 70% UT for 25/30 cycles at 0°.	The quality of the power supply must be the same as in a home, business or hospital environment.
IEC 61000-4-39	8 A/m (30 kHz, CW) 65 A/m (134.2 kHz, PM 2.1 kHz) 7.5 A/m (7.5 kHz, PM, 50 kHz)	8 A/m (30 kHz, CW) 65 A/m (134.2 kHz, PM 2.1 kHz) 7.5 A/m (7.5 kHz, PM, 50 kHz)	The magnetic field at the network frequency must be the same as for a home, business or hospital environment.

Test specifications for port immunity from box to RF wireless communication equipment:

TEST FREQUENCY (MHz)	MODULATION	IMMUNITY TEST LEVEL (V/m)
385	18Hz pulse modulation	27
450	FM, ±5 kHz deviation 1kHz sine	28
710	Pulse modulation 217 Hz	9
745		
780		
810	Pulse modulation 18 Hz	28
870		
930		
1720	Pulse modulation 217 Hz	28
1845		
1970		
2450	217 Hz pulse modulation	28
5240	Pulse modulation 217 Hz	9
5500		
5785		


XII - BEFORE CALLING DRIVE DEVILBISS HEALTHCARE FRANCE CUSTOMER CARE SERVICE

PROBLEM	INDICATOR / SYMPTOM	POSSIBLE CAUSE	SOLUTION
The patient lift is not stable		The mast-base junction has some play in it	Refer to the installation of the mast on the base in the Installation section
Loud noise during a movement, significant friction	Noise coming from the wheels, difficulty pushing the patient lift	Dust in the wheels	Remove and clean the wheels
	The noise comes from the pivots	Dust in the pivots	Remove and clean the joint
The screen is lit, a button is pressed but the cylinder does not lift	Rotation of the electric motor but no lifting of the cylinder rod	The free rod system was engaged after hitting an obstacle	Release the patient lift from the obstacle and guide the rod back to its normal position in the cylinder body.
The screen is black	The screen does not light up	Emergency stop button engaged	Unlock the emergency stop button
		Control box in energy-saving mode	Press any button to exit energy-saving mode
		Discharged battery	Charge the battery (direct or indirect charge)
		The remote control is incorrectly connected	Press a button on the control box to exit the energy-saving mode. Check the remote control connection in the control box
		Faulty control box	Contact the Drive DeVilbiss Healthcare France customer care service
Lifting takes place in down-graded mode	The cylinder lifts but at a reduced speed	Battery has a low charge	Charge the battery
Lifting function does not work	Rotation of the electric motor but no lifting of the cylinder rod	The free rod system was engaged after hitting an obstacle	Release the patient lift from the obstacle and guide the rod back to its normal position in the cylinder body.
		Unusual cylinder noise	Contact the Drive DeVilbiss Healthcare France customer care service
	No motor rotation	Protection against overheating has been activated	Disconnect and allow the system to cool for 20-30 minutes
		The overload protection has been activated. The signal is displayed on the screen and on the remote control	Reduce the load to be lifted
		The connection cable between the control box and the cylinder is incorrectly connected	Check the cable connection

PROBLEM	INDICATOR / SYMPTOM	POSSIBLE CAUSE	SOLUTION
The battery will not charge at the charging station	When plugged into the wall socket, the charging station indicator does not display green without a battery.	The transformer cable is incorrectly connected to the charging station	Check the cable connection
		The charging station is faulty	Contact the Drive DeVilbiss Healthcare France customer care service
	The charging station displays a green light without a battery, but does not change to orange when the battery is connected	The battery is incorrectly connected to the charging station	Check the battery connection
		The battery is defective	Contact the Drive DeVilbiss Healthcare France customer care service
	The charging station is green after the charging time, but the battery is discharged once it is placed in the control box.	The battery can no longer be charged	Contact the Drive DeVilbiss Healthcare France customer care service

 **Contact the Drive DeVilbiss Healthcare customer care service if the above solutions do not work.**

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