

User Instructions

Standard Floor Fixed Folding Support Rail

M11260 Long 760mm (30+)
M20108 Short 550mm (21½+)



Maximum User Weight

General Safety Warning

This product should only be fitted by a suitable qualified person.

Please ensure these instructions are fully read and understood before installation/ use. Retain these instructions in a safe place for future reference.

Intended Use

The Floor Fixed Standard Support Rail is a high quality steel rail which provides a stable support to assist the user raise and lower themselves for toileting. Suitable where the wall is unsuitable for a wall fixed version. The rail can be raised back when not in use or when transferring from the side.

Safety Precautions

This device should only be secured to sound wooden floors on timber joists or concrete floors. It should only be installed by a suitably qualified person.

Periodically check the stability and rigidity of the product and, where applicable, that screws are tightened and in good condition.

Do not exceed maximum user weight.

Installation Instructions

1. Place stand in required position. Ensure that the rail can be raised and lowered without interference from other fixtures and furniture and that transfer can be carried out safely.

2a. For floor boards . make sure screws go through floorboards into joists. Drill 5mm pilot hole in each of the 4 screw positions. Screw down using appropriate fixings.

2b For Concrete Floors . using a masonry drill, spot through the 4 holes in the base plate. Remove stand and drill 4 holes as appropriate for the fixings. Remove dust, insert plugs and replace stand. Secure using suitable fixings.



3. Using the black hand-wheel on the side of the column, release the support rail and slide up/down the column. Re-clamp the rail securely at the appropriate height using the black hand-wheel.

4. When fitted, check the product is securely fixed to the wall and that the arm can be raised and lowered without interference. Also check that the friction lock is effective and that the support rail will hold itself at any angle during movement (i.e. it does not drop down on it's own).

Cleaning Instructions

Clean all parts regularly, wiping down with mild detergent /disinfectant solution.

Manufactured for

NRS Healthcare

Sherwood House, Cartwright Way,
Forest Business Park, Bardon Hill,
Coalville, Leicestershire.
LE67 1UB UK

Tel. Main Office: 0330 100 00 98

Tel. Customer Services: 0845 120 4522

Email: customerservice@nrs-uk.co.uk.

This product is from a range of over 3000 daily living aids available from NRS Healthcare. To view the full range visit www.nrs-uk.co.uk or request a catalogue from the above number.