



User Instructions

NUVO™ Slatted Bath Seat

Nuvo™ Slatted Bath Seat Range

L71055 / L71079	Bath seat 6/8" Combo
L71067	Bath seat 10/12" Combo
L71092	6/8" Leg inserts
L71109	10/12" Leg inserts + base
L71110	Spacers
L71122	Seat Top

Maximum
User
Weight



Please ensure these instructions are fully read and understood before use. Retain in a safe place for future reference.

Intended Use

The Nuvo™ Slatted Bath Seat is designed to provide a safe seat for use when bathing. It can be used on its own or in conjunction with Nuvo™ Bath Boards to assist the user in and out of the bath.

Prior to fitting, please ensure:


- The user does not exceed the 190kg (30 stones) weight limit
- The bath surface is clean, dry and free from oil or soap residue
- The correct size of bath seat has been selected for the client
- Do not use in a highly sculptured or unusual shaped bath, as the seat may not fit properly.
- Do not use in a bath with a textured base
- Do not stand on the seat
- Always ensure that the sucker feet are secure before bathing commences

Installation Instructions

1) Screw leg inserts onto base in a clockwise direction. If an 8+seat is required, screw spacers onto base first, as shown . see fig. 1

If a 12" seat is required screw in spacers place base frame over spacers, screw in leg inserts then pull base frame firmly down to locate over each rounded foot

Fig. 1



2) Place the Bath Seat in the desired position on the flat base of the bath and push down firmly to ensure that the sucker feet have adhered to the bath surface.

3.) Ensure that the slats are positioned parallel to the width of the bath and that all sucker feet are secure.

Removing

Before lifting, release the tabs on the sucker feet.

Care of the Nuvo™ Bath Seat

The Nuvo™ Bath Seat is made from corrosion free parts. It can be cleaned with non-abrasive cream cleaners or up to 80 C in industrial cleaners. Please complete a visual inspection of the seat every six months, checking for wear and tear.

Safe Use of the Nuvo™ Bath Seat

6”(150mm) / 8”(200mm) versions

It is recommended that the 6+and 8+Nuvo™ Bath Seats are used in conjunction with a Nuvo™ Bath Board, for optimum safety. The recommended safe method of using a Bath Board and Seat together is as follows:

1. Start with your feet flat on the floor and your back to the bath
2. Step backwards until you can feel the edge of the bath against the back of your legs
3. Sit on the board and move your bottom back towards the middle of the bath board
4. Lift your legs one at a time over the bath edge, into the bath
5. Move your bottom across the board until you are seated in a central position. You may need to use your arms to lift your bottom and move sideways in small stages

You are now ready to use your arms to lower yourself down on to the Nuvo™ Bath Seat
Reverse the procedure for getting out again
the top of the bath seat is level with the top of the bath.
In other cases, the user should step into the bath (where appropriate) and sit down onto the seat, possibly using a grab rail for stability.

10”(250mm) / 12”(300mm) versions

1. It is recommended that the safest way to use the 10/12+Nuvo™ Bath Seat is to use the same method as for the Nuvo™ Bath Board above (stages 1-4).
2. Reverse the procedure for getting out again.

Please note, this method is only recommended in situations where the top of the bath seat is level with the top of the bath. In other cases, the user should step into the bath (where appropriate) and sit down onto the seat, possibly using a grab rail.

Still Having Difficulty?

Here are some ideas that may help you:

Problem	Suggestion
Difficulty swinging legs over bath edge	Try a manual leg-lifter (NRS Code L33194)
Difficulty getting into position on the board or seat	Try a grab rail, mounted on the wall above the bath. NRS have many different sizes and styles. Call Customer Services for details.
Difficulty reaching down to wash lower legs and feet or up to wash hair or scalp	NRS have a range of long-handled personal grooming aids to assist with this. Call Customer Services for details.

IMPORTANT

For information on how to obtain the User Instructions in other languages and formats for use by people with visual, reading or cognitive impairments, please contact: info@nrs-uk.co.uk or telephone 0845 120 4522.

Additional copies of this user instruction are available on our website www.nrs-uk.co.uk.

Replacement Products & Spare Parts

Please contact NRS Customer Services.

Tel: 0845 120 4522

Fax: 0845 121 8112

Manufactured in China exclusively for

NRS Healthcare

Sherwood House,
Cartwright Way,
Forest Business Park,
Bardon Hill,
Coalville,
Leicestershire.
LE67 1UB UK

Tel. Main Office: 0330 100 00 98

Tel. Customer Services: 0845 120 4522

Email: customerservice@nrs-uk.co.uk

www.nrs-uk.co.uk 