

Usage

- · Be sure to read all information and directions before use.
- The Quick Rail is designed to assist your natural motion, as when stepping into or out of a shower.
- The Quick Rail is not designed to hold or support your full body weight.
- Proper placement of the Quick Rail is according to individual use.
- · Please be aware of your own ability and limitations in strength and range of motion.
- The *Quick Rail* can be used where drilling to mount bars is not possible; for example in fiberglass bath enclosures, or glass shower panels.
- · Do not attach to moving surfaces.
- The Quick Rail is suitable for rental situations, or to carry with you for travel.

Attaching the Quick Rail

- Quick Rail is designed for use only on smooth and non-porous surfaces; tile, glass, porcelain or fiberglass. Not for use on wall board, wood, painted or papered surfaces.
- Check the condition of the wall or tub surface before positioning the *Quick Rail*. Support surface should be free of defects and be in sound condition.
- · Clean the surface thoroughly with alcohol and allow it to dry before attaching the Quick Rail.
- · Never mount in a position when the Quick Rail will be in direct sunlight.
- · Check the condition of the suction cups. Do not use if damaged, cracked or split.
- With both blue levers raised, place the suction cups so that they lie firmly and completely on a smooth surface. Not for use on mosaic style tiles.



- Do not place over grout lines as illustrated at left.
- Press on one lever on one end, then the other. You should feel resistance as it tightens.
- Check attachment before each use.
- To remove, lift each of the 2 levers. If the suction remains, lift one edge with a finger. Do not use a sharp edge, as damage to the suction cup is possible.
- Clean Quick Rail regularly with non-abrasive alcohol based bathroom cleaner.

Liability warning

Quick Rail is designed for use only on smooth and non-porous surfaces. Always check attachment before use. Before installing, check the stability of the wall. All tiles must be firmly attached to a solid backwall. One piece molded inserts must be firmly adhered to the wall. Quick Rail is sold under the condition that the user assumes all risk of injury due to use.



Days Healthcare GmbH & Co. KG Oberbecksener Strasse 68 32547 Bad Oeynhausen

(t) + 049 (0)5731.78 65 - 0

(e) info@dayshealthcare.de

(f) + 049 (0)5731.78 65 - 20

(w) www.dayshealthcare.de