

User Instructions for Childs Parallel Bars

Maximum User Weight: 63kg

THIS PRODUCT SHOULD ONLY BE USED UNDER THE SUPERVISION OF A TRAINED THERAPIST

Please read these instructions fully before using this product and retain for future reference.



These bars can be used fixed to the floor, or free standing. When fixed to the floor use 8mm fixings installed by a qualified engineer.

NOTE: It is recommended that these Parallel Bars are fixed to the floor when in use. Extra vigilance is required when they are used when not fixed to the floor.

When assembling or moving the assembled bars, always use 2 persons and make sure the locking pins are fitted in place properly.

When folding or disassembling the bars, always use 2 persons and remove the locking pins carefully noting the possibility that the bars may drop suddenly with the pins removed.

Height Adjustment

Insert the 4 clips at required height. Remove the clips, lower hand rails to fold away.

Care and Cleaning

No sterilisation is required but bars should be wiped over between users with a mild detergent solution or surface disinfectant to avoid cross contamination.

Cautions

When used free standing be careful not to allow the user to push on the bars in the direction of travel as the bars may slide forwards. Do not allow the user to push on the bars in a sideways direction as the bars may tip over.

Do not use or store the bars in a damp environment.

Have the bars checked annually by a competent engineer.

Manufactured by G Shaw & Sons Ltd, Nottingham

