




Support Rail with Concealed fixings

Fitting and user instructions

Congratulations on purchasing your support rail.

The Support Rail offers a safe hand hold providing:

- a support to push or pull against when standing up or sitting down.
- a steadying support & balance when standing or dressing.
- a firm grip whilst transferring from one position to another.
- aesthetically pleasing concealed fixings.




Fitting

WARNING: Before fitting, ensure that there are no concealed electrical cables or pipes where the fixing holes are to be drilled. If installing a metallic rail in a bath/shower area the rail may need to be earthed (IEE wiring regulations); not all rails are suitable for use in a wet area.

1. Ensure the wall and fixings being used are suitable. Nymas are unable to specify fixings due to the number of fixing/wall type combinations.
2. Position the Support Rail on the wall either horizontally or at an angle to suit and mark the hole locations.
3. Drill the holes to the correct depth for the wall and fixing type used.
4. Ensure the rail is securely fastened to the wall and fit flange covers to conceal the fixings.

If in doubt about anything on this sheet consult a suitably qualified professional.



Using the Support Rail

- Before first use, check that the Support Rail is securely fitted to the wall.
- Use the Support Rail to assist; **NEVER** apply your full weight to the support rail.
- Regularly clean the Support Rail with a non-abrasive proprietary cleaner.
- All fixings and bolts should be checked on a regular basis by a competent person.

