

Fitting and Operational Instructions – continued;

6 - Attach the shoulder straps to the hoist spreader bar. (The middle loop on each shoulder strap will give a reasonably upright sitting position.) A more reclined position can be achieved by using a longer loop. For a more upright posture choose a loop closer to the sling.

7 – Attach the head support panel straps to the spreader bar, this must be done after the shoulder straps have been fitted. To achieve a comfortable position in which the client's head is correctly supported the straps can be adjusted by opening the cam buckle on the strap, pulling the strap through until the correct length is achieved; the cam buckle should then be closed securely. Note: Inspect the cam buckles, regularly, for signs of wear or damage.

8 - Attach the leg straps to the hoist spreader bar. The middle loop on each leg strap will give a reasonably natural sitting position. The client's knees can be raised by selecting a loop closer to the Sling.

9 - Ensure that all loops are fully attached to the spreader bar, and that the sling is well positioned and not attached in any way to the seat or other obstacle.

10 - Ensure client is secure and comfortable.

11 - Before commencing the hoist, carers must be in the pool to ensure that the transfer into the water is achieved safely and that all normal safety considerations are observed.

12 - Raise the hoist spreader bar to lift the client just off the surface. Check that all loops are connected and the client is secure and comfortable prior to completing the transfer.

WARNING: Do not attempt to adjust or detach any strap while the client is being hoisted.

Carers must support the client whilst straps are detached. Note the Head support straps must be detached before the shoulder straps.

Fitting from the Lying position

Roll the client onto the sling, ensuring that the Bottom edge of the Head Support panel is level with the base of the neck, and the bottom edge is positioned just below the coccyx. Follow instructions from point 3 of the seated position instructions.

NRS Healthcare

Sherwood House, Cartwright Way,
Forest Business Park, Bardon Hill,
Coalville, Leicestershire.
LE67 1UB, UK

Tel. Main Office: 0330 100 00 98
Tel. Customer Services: 0845 120 4522
Email: customerservice@nrs-uk.co.uk

www.nrs-uk.co.uk



User Instructions

Deep Back Support Sling Range

Polyester Mesh

Sling Guide					
Size/Colour Code	Ex. Small	Small	Medium	Med/large	Large
Sling Code	M07232	M07244	M07256	M07268	M07279
Client Measurements Head to Sacrum	500 - 600 mm 19½"- 23½"	600 - 850 mm 23½"- 33½"	850 - 1050 mm 33½"-41"	930 - 1130 mm 36½"- 44½"	1050 - 1200 mm 41"-47"

Maximum User Weight



General Safety Warning

Please ensure these instructions are fully read and implemented. Failure to do so may result in injury to the user. Retain in a safe place for future reference.

In the interest of safety; Persons using this equipment should have received appropriate training or instruction. Hoisting transfers should only be attempted following a comprehensive risk assessment and should be performed in accordance with local Moving & Handling Guidelines.

Intended Use

The Deep Back Support Sling range is designed to suit a wide range of users and transfer needs, including bathing, and use in swimming or hydrotherapy pools. They can be used on hoists fitted with a conventional spreader bar (i.e. with hooks at either end from which looped straps may be suspended)

If you are unsure about compatibility with the spreader bar contact NRS customer services, contact details overleaf.

Cleaning Instructions:

Note: The head support stiffening bars (Positioned in the head support section of the sling) **MUST** be removed before the sling is washed.

- Machine wash warm at up to 70°C
- For effective decontamination, wash at greater than 65°C for at least 10 minutes
- Do not use bleach, biological detergents or detergents containing bleach.
- Do not spot clean with strong chemicals
- Do not iron.
- Do not dry clean
- Tumble dry only on 'cool' setting or air dry at very low temperatures.

Special Note:

After use in a swimming or hydrotherapy pool slings must immediately be rinsed thoroughly in cool, tap water, and dried.

Inspection Information

The following notes are a guide to the recommended inspection procedures for slings.

If carried out thoroughly, at regular intervals, lifting and hoisting operations can be carried out with confidence, and the risk of injury due to sling failure is significantly reduced.

Before every lift, visually inspect the slings structural integrity. Use the diagrams as a guide. If there are any signs of damaged fabric or stitching do not use.

If you are unsure, do not use the sling, withdraw from service immediately. The product should then be either discarded, or clearly labelled:



“Not for use” and assessed at the earliest opportunity by a competent individual.

Sling Inspection procedure;

The following procedure is recommended to help ensure that lifting and hoisting operations can be carried out with confidence and the risk of accident due to sling failure is minimised.

- Open the sling out fully on a flat surface.
- Check all retaining straps for cuts, nicks or degradation of any kind.
- Carefully check all stitching.
- Check the material near to any attachment points.
- Pay particular attention to the areas of high stress as shown in diagram.
- Checks should be carried out to both sides of the sling in the same manner.
- Discolouration may also indicate fabric degradation.


CAREFULLY EXAMINE ALL STITCHING




CHECK:

- All Stitching **MUST** be intact.
- There should be no Material Degradation around stitching.

STRAPPING INSPECTION



A visible nick – even of only 1mm **SHOULD BE DEEMED UNSAFE**



A few loose strands are **ACCEPTABLE**

Fitting and Operational Instructions

Before every lift, visually inspect the sling to ensure it is fit for purpose. Refer to the Inspection Routine section of this guide for the recommended method. If in doubt do not use the sling.

In the interest of safety; Persons using this equipment should have received appropriate training or instruction. Hoisting transfers should only be attempted following a comprehensive risk assessment and should be performed in accordance with local Moving & Handling Guidelines.

Fitting from the seated position

1 - Place the Sling behind the client with the bottom edge of the head support panel level with the base of the client's neck.

2 - Position the sling by pulling the bottom edge briskly from left to right and back again, the bottom edge of the Sling can be positioned 25-50 mm under the client's buttocks.

3 - Pull the leg straps forward parallel with the client's upper leg. (It may be helpful to “block” the client's knee with one hand while this is being done).

4 - Pull each leg strap under the adjacent upper leg, so that both straps emerge between the client's legs.

5 - Each leg strap should now be positioned comfortably under the client's leg. If not, go back to step 3 and repeat.