

School Chair Footrest



Bath Institute of Medical Engineering

Wolfson Centre

Royal United Hospital

Bath BA1 3NG

Tel: 01225 824103 Fax: 01225 824111

Email: info@bime.org.uk

Website: www.bime.org.uk

Registered charity 256335

School Chair Footrest

The footrest is designed to clamp to the legs of standard school chairs. To assemble and fit to your chair you will need two 13mm spanners. (Not supplied.)



Assembly of parts of school chair footrest



Brackets:

20mm diameter

Determine height of footrest required:-

Adjust the legs of the footrest to the correct height for the child. The height can be adjusted between 115 and 160mm by depressing the spring loaded button. Slide the leg until the button inserts and locks in the hole required for the height selected.

Ensure the button has popped back fully through the hole, locking the leg at the selected height.



Button to adjust height of legs of footrest


To assemble and attach the footrest:-

✎ Slide the brackets into the footrest.

✎ You will now need to unscrew the fixing assembly, in order to slide the bracket onto the chair legs. You may need to ease open the brackets using tools to allow the fixing to be fitted to the legs. Before securing the fixing assembly ensure the footrest is at the correct height and that the platform is level and stable.

✎ Cover the ends of the fixing assembly with the plastic sleeves provided.

See order of fixing assembly:-



Tightening the fixing assembly

✎ Before use check that the footrest is securely fitted with a level platform and with all legs fully resting on the floor. Ensure that the footrest platform can be folded away under the chair.



Footrest Positions:-

In use the footrest should be resting squarely on the floor. When not in use or for transfers on to the school chair the footrest can be folded under the chair and is supported by the fixing assembly. Always fold by holding the front of the footrest.



Please note: *Our footrest will fit infant size chairs as shown. If a footrest is needed for a junior size classroom chair, please contact us for details of larger size clamps and legs to use with your school chair footrest.*

CAUTIONS:

 **Ensure the footrest is stable before use.**

 **Always ensure the height adjustment button is protruding well through the hole when adjusted and during use**

 **Check fittings periodically to ensure footrest is safely secured to school chair**

 **Be aware when folding or touching the footrest that you do not run the risk of pinching or trapping of hands / fingers**

 **Do not stand on footrest using it as a footstool**

WARRANTY

The product is guaranteed from date of purchase for one year.