

Technical Information:

Cushion/ Quality	Sit Right in Bed
Cover	Simulated Lambswool
Foam Inner	Combustion modified polyurethane foam
Dimensions (Max)	59 x 53 x 10cm

Product Contents:

1 x Sit Right in Bed

Contact Details: If you require any further information please contact us at:

Thorpe Mill Ltd
Unit 1a Aireside Business Park
Royd Ings Avenue
Keighley
BD21 4BZ

Tel: 01535 682630
Fax: 01535 682639

sales@thorpemill.co.uk

Please take time to read the contents of this information guide. Any actions that are inconsistent with this guide, are done so at the operators risk and Thorpe Mill Ltd will not be liable for injury or damage. In addition to the information in this guide, operators should use the product specified with due consideration of applicable international, national and local regulations, guidance and policies. Safe and correct use of this product is solely at the discretion of the operators.

If a serious incident does occur during use, this should be reported to Thorpe Mill Ltd using the above contact details.

Instruction For Use Sit Right in Bed (BRA)



Product Name: Sit Right in Bed (BRA)

Intended Purpose and Design: The Sit Right in Bed is a support wedge, designed to make sitting in bed more comfortable. With its armchair style back and arm rests, the Sit Right in Bed helps the user stabilise and sit up in bed with ease. Comfortable and easy to use, the Sit Right in Bed is made of combustion modified polyurethane foam and features a simulated lambswool cover, with side pocket.

Operator Qualifications: Those responsible for the storage, transportation, usage and cleaning of the Sit Right in Bed should ensure they read these instructions in full. The Sit Right in Bed should only be used for its intended purpose, as specified in this instruction guide.

If the Sit Right in Bed is to be used in a care/medical setting, use of the Sit Right In Bed should be determined by a qualified health professional, following a comprehensive and holistic assessment using both formal and informal assessment methods. In addition to the information in this guide, operators must adhere to any local manual handling and health and safety policies.

Instructions for Use (IFU):

Prior To Use: Upon receipt of the Sit Right in Bed all packaging should be removed and visual checks should be carried out to ensure:

- The contents are complete
- The Sit Right in Bed has not been damaged in transit

Storage: When the Sit Right in Bed is not in use it should be stored safely to avoid damage. If in a care/medical setting, the Sit Right in Bed should be stored in a secure location away from the general public. To avoid damage to the Sit Right in Bed, never store other items on top of it, do not store it directly on the floor and do not store it next to radiators or other heat sources.

During Use:

- Place the Sit Right in Bed on top of the mattress (and any relevant bedding), ensuring it is positioned safely, in the middle of the bed against the head rest
- If the user has difficulties with mobility/stability, ensure that the bed has appropriate bed rails in place

Aftercare: The Sit Right in Bed is a reusable product, therefore, the following processes should be undertaken to allow its re-use.

Visual Inspection: After use, the Sit Right in Bed should be visually checked to see whether it is in good condition and for signs of material degradation or damage. Check for any signs of tearing and/or punctures, check all the seams for any signs of splitting, check the zip(s) for any signs of damage and check for any signs of permanent contamination.

Unzip the Sit Right in Bed and check the foam inner for any signs of contamination. If following the visual inspection, there are any signs of contamination, the Sit Right in Bed should be withdrawn from use until the cover and/or the inner foam is replaced. If the cover or foam inner is replaced, then it is important that the correct corresponding cover/foam inner is re-ordered, as Sit Right in Bed covers must not be used on other foam inners.

Contaminated covers/inner foams, should be disposed of safely and in line with local infection prevention and control policies, if applicable. If the Sit Right in Bed is in good condition, then it should be cleaned as per the instructions below and stored in line with the storage instructions.

Cleaning: The Sit Right in Bed cover is machine washable up to 30°C. The Sit Right in Bed can be flat dried in an indoor, clean environment. Do not tumble dry and do not bleach. Sit Right in Bed covers must be thoroughly dried before re-use. Always ensure the correct cover is returned to its corresponding foam inner. Covers must not be used on other inner foam inserts.



General Warnings Precautions and Limitations:

- Only use the Sit Right in Bed for its intended purpose
- Take care not to rip/puncture the cover, to avoid contamination
- Do not place sharp objects such as scissors, needles, syringes, scalpels or other sharps on the Sit Right in Bed
- Advise operators, users and visitors to remove or cover hand jewellery including rings with sharp edges
- If applicable, take extra care when using medical equipment near the Sit Right in Bed
- Do not drag along the floor, do not transport in roll cages unless completely protected from sharp edges and do not scrape against walls, door frames, door catches or locks
- Avoid any protruding and sharp objects
- The cover must be thoroughly washed and dried between and after cleaning. The cover is more susceptible to physical damage when wet and will remain susceptible for a period of time after being dried
- Ensure the Sit Right in Bed is placed on top of the mattress (and any relevant bedding), ensuring it is positioned safely, in the middle of the bed against the head rest
- If the user has difficulties with mobility/stability, ensure that the bed has appropriate bed rails in place
- If the Sit Right in Bed is to be used in a care/medical setting, use of the Sit Right in Bed should be determined by a qualified health professional, following a comprehensive and holistic assessment using both formal and informal assessment methods

