



The **MPPL-FMAT** is a simple to use paging system comprising of a transmitter, a tone and vibrate pager and a large non slip pressure mat.

## SETUP:

**The battery is already inside the transmitter, but not connected. To begin slide down the bottom half of the case and connect the battery, observing polarity.**

The pager requires 2 x AA batteries (already in place).

You can now turn the pager on by selecting the alarm type and pushing the side switch to the 'on' position. 'Vib only' is a vibration setting, 'Vib + vol' is a low volume and vibrate setting and 'Vib + Vol' is a loud tone and vibrate setting. On any selection the pager needs to be switched to off to reset.

## Pressure Mat

The pressure mat has a large socket and the transmitter has a long lead with a plug on the end. Insert the plug into the socket they should click together. Position the mat on the floor appropriately. Slide the power switch on the side of the transmitter and also the pager. Standing on the mat will activate the pager.



Medpage T/A EasylinkUK

3 Melbourne House

Corby Gate Business Park

Priors Haw Road

Corby

Northants

NN17 5JG

Tel 01536 264869