

When you are in south or east area of UK, the pill dispenser would receive RCC signal of Germany instead of UK. In such case the “**DCF**” symbol will be shown on the screen, as below:



The clock will be 1 hour ahead of actual UK time.
Below procedures show how to adjust this time difference.

Unlock the pill dispenser and open the top cover.
Press Button 2 first, hold and press Button 1 simultaneously for a few seconds until “H 00” is shown on the screen.




Press Button 3 once to adjust the time difference to **-01** hour,
“H - 01” will be shown on the screen:
(In case you press Button 3 for more than once, you may press Button 2 to
go backwards)




Press Button 1 once to quit time difference adjustment,
The original clock time will be shown on the screen:



Press Button 3 and hold until the RCC symbol  is flashing on the screen, which means the pill dispenser is scanning for RCC time signal again.



Once the RCC signal is received, the RCC symbol  will stop flashing and the correct time (with adjustment of -1 hour) will be shown on the screen:

