

User Instruction for the Mobile Cream Split Level Top Table B45



1. Introduction

The Split Level Top Table is designed to aid the patient, by being suitable for use over a bed or chair – offering a high quality, versatile service. The top is hard wearing and easy to clean, and the retaining table edges prevent items from slipping off. There is a tilt-facility that is simple to operate and the product is adjustable in height. It features castors for easy movement.

Product Dimensions	
Large table top size	610mm x 458mm
Small table top size	458mm x 150mm
Table top height (max)	940mm
Table top height (min)	630mm
Footprint	685mm x 400mm

2. Contents of the Box

- Split Level Top Table (B45)
- Extended Legs & Handwheels
- Castors
- User Instructions (this leaflet)

3. Assembly Instructions

- Carefully unpack and examine the product to ensure it is complete.
- Before use check for any signs of damage e.g. cracked weld, flaking of paint, bent tube. If any signs occur please contact your supplier to exchange and refrain from using the product. **DO NOT ATTEMPT TO USE IF DAMAGED.**

4. Correct Use

- The table is designed to help the user by being used as a dining, reading, or writing aid.
- The table is not intended to support full body weight so refrain from attempting to lift full body weight via the product
- The table should be placed either against a bed or chair.
- Do not attempt to overload the table.

Adjustable Height:

- Unscrew both hand wheels.
- Raise/lower to required height.
- Retighten the hand wheels.
- Tilt table top to angle required by raising or lowering on each side of the table.

User Instruction for the Mobile Cream Split Level Top

Table

B45

- Tighten both hand wheels to lock top into desired position.

5. Cleaning

- Use mild detergent with warm water to wipe down framework on at least a weekly basis more frequently if required.
- Do not use an abrasive cleaner as this will damage the frame and nullify your warranty.
- **The beech table top is not suitable for machine cleaning and should be removed and cleaned as above if the frame is to be cleaned by mechanical means.**
Use of power/jet washing & industrial washing machines:
- If either of the above is used to clean equipment care must be taken since high-pressure water can be penetrative.
- Care should be taken that frames are not placed too close to high-pressure water jets.
- After cleaning ensure any residual water is removed from the framework as it could cause permanent rusting of the frame (leave frames open in a warm environment for several hours to ensure all internal moisture evaporates).
- We cannot offer advice on industrial washing machines with regards to specific machine settings, as these will be dependant on the manufacturer's instructions.

6. Maintenance Requirements

- The product should be checked on a regular basis for damage or wear and tear – this includes flaking or cracking paint. Any undue signs of damage should be noted and reported to the table supplier in the first instance. Should damage be evident, then the user should refrain from using the product until advice has been sought from the supplier..

7. Warning

- This product is intended for use as a dining, reading or writing aid – it does not support heavy weights.
- We recommend a full risk assessment should be carried out before use.

8. Warranty

This product carries a 12 month warranty from the date of purchase against manufacturing defects. Misuse or damage is not covered. This warranty does not affect your statutory rights. As manufacturers we cannot specify expected product life as this is based on the level of usage.

9. Disposal Instructions

- Appropriate hygiene standards should be applied before disposal to reduce the risk of infection.
- Materials used are essentially non hazardous and can be disposed of appropriately.

10. Contact Details

We have attempted to cover most eventualities, however if you have any specific queries, please contact:

Cefndy Healthcare
Cefndy Road
Rhyl
Denbighshire
LL18 2HG
Tel: 0044 1745 343877 (Quality Assurance Dept)
E-mail: cefndy.sales@denbighshire.gov.uk