

## User instructions & technical data sheet Standard Sero Cushion Bonyparts Cut Out (SPCBY)



Maximum user weight: 100kg

Risk category: Low (Please Note: Indication of pressure ulcer risk is not a prescriptive evaluation. Users need to be individually assessed by a tissue viability expert).

User instructions: Sit with profile “egg box” facing upwards and cut-out towards back of seated area. Do not wash foam or get it wet. Cover can be removed and hand washed.

Polyurethane, high-density foam with zipped, removable cover. Meets BS5852.

Manufactured by:


Putnam Health Co Ltd, Eastern wood Road, Langage Industrial Estate,  
Plympton, Devon, England, PL7 5ET

Operations Manager:

Hedley Putnam

## Independent pressure mapping of Standard Sero Cushion Bonyparts Cut Out (SPCBY) V's control cushion.

Standard Sero Cushion Bonyparts Cut Out (SPCBY)



Control cushion

