



P42434: NRS EasyFit Bedguard - Single

P78222: NRS EasyFit Bedguard - Double

The NRS EasyFit Bedguards are designed to prevent users from rolling out of bed.



The NRS EasyFit Bedguard is designed to prevent users from rolling out of the bed.

The Bedguard can be used on divan and slatted beds, from single to king size.

It is designed to be very easy to install to a bed, without needing securing straps.

It has a telescopic side rail, which can be extended or retracted for easy access on and off the bed.

The side rail can be removed to facilitate making the bed.



P42434: Single sided model, to protect users from rolling out of one side of the bed only.




P78222: New double-sided model for users who are at risk of rolling out of either side of the bed.

Product Features





Assembled in 4 easy steps, with no tools or straps required. Features an adjustable base to fit different bed widths.



Telescopic side rail for easy access on and off the bed. Mesh cover provides additional safety, and includes pockets for storage of books etc.

Product Dimensions



For more information or to place an order, please contact your local Regional Sales Manager or the NRS Healthcare Customer Service Team:

Telephone: 0345 121 8111

Email: customerservice@nrshealthcare.co.uk

www.healthcarepro.co.uk or www.completecareshop.co.uk

For international enquiries, please telephone: +44 (0)1530 232 292

or email: exportsales@nrshealthcare.co.uk

