

General Information

Please ensure these instructions are read and understood before use. Retain them in a safe place for future reference.

Intended Use

This one-piece lightweight doorframe ramp provides wheelchair users with a convenient solution for door thresholds. Clearance level of up to 70mm (2¾"). Light to carry and can be easily stored when not in use. Made from aluminium with a durable slip resistant track for grip and safety. The ramp is for indoor use.

To view a video showing how to install and use this product go to:

<http://www.nrs-uk.co.uk/doorframe-ramp.html>

Dimensions:

Length: 900mm (35")
Width: 690mm (27")
Floor Clearance: 70mm (2¾")
Max Weight Load: 130kg (21st)

Cleaning Instructions

The Doorframe ramp can be cleaned with domestic cleaning materials without abrasive or corrosive properties. Solvents and chlorine should not be used.



Product Guarantee

NRS Healthcare, guarantee this product for a period of 12 months from the date of purchase against faulty materials and workmanship. This guarantee does not cover, abuse, misuse, interference or tampering with any part of the product. Any guarantee issued by NRS Healthcare does not affect your statutory rights.

Manufactured in the UK and Distributed by:

NRS Healthcare

Sherwood House, Cartwright Way
Forest Business Park, Bardon Hill
Coalville, Leicestershire.
LE67 1UB UK

Tel. Head Office: 0330 100 0098
Tel. Customer Services: 0345 121 8111
Email:
customerservice@nrshealthcare.co.uk

NRS Healthcare is a Leading Supplier of Community Equipment and Daily Living Aids, to view our full range, visit our website: www.nrshealthcare.co.uk

General Information

Please ensure these instructions are read and understood before use. Retain them in a safe place for future reference.

Intended Use

This one-piece lightweight doorframe ramp provides wheelchair users with a convenient solution for door thresholds. Clearance level of up to 70mm (2¾"). Light to carry and can be easily stored when not in use. Made from aluminium with a durable slip resistant track for grip and safety. The ramp is for indoor use.

To view a video showing how to install and use this product go to:

<http://www.nrs-uk.co.uk/doorframe-ramp.html>

Dimensions:

Length: 900mm (35")
Width: 690mm (27")
Floor Clearance: 70mm (2¾")
Max Weight Load: 130kg (21st)

Cleaning Instructions

The Doorframe ramp can be cleaned with domestic cleaning materials without abrasive or corrosive properties. Solvents and chlorine should not be used.



Product Guarantee

NRS Healthcare, guarantee this product for a period of 12 months from the date of purchase against faulty materials and workmanship. This guarantee does not cover, abuse, misuse, interference or tampering with any part of the product. Any guarantee issued by NRS Healthcare does not affect your statutory rights.

Manufactured in the UK and Distributed by:

NRS Healthcare

Sherwood House, Cartwright Way
Forest Business Park, Bardon Hill
Coalville, Leicestershire.
LE67 1UB UK

Tel. Head Office: 0330 100 0098
Tel. Customer Services: 0345 121 8111
Email:
customerservice@nrshealthcare.co.uk

NRS Healthcare is a Leading Supplier of Community Equipment and Daily Living Aids, to view our full range, visit our website: www.nrshealthcare.co.uk