

Fitting instructions for Evacuation Sledge (EVSG).

This Evacuation Sledge is designed to be stored in its bag on the wall for easy access in case of emergency

- Pull pocket strap to release Sledge.
- Lay Evac Sledge fully extended on the floor cotton side up
- Lay patient on the Evac Sledge.
- Ensure patients feet are in the foot pocket.
- Tighten security belts over patient.
- Pull the patient to safety.

Fitting instructions for Evacuation Sledge (EVSG).

This Evacuation Sledge is designed to be stored in its bag on the wall for easy access in case of emergency

- Pull pocket strap to release Sledge.
- Lay Evac Sledge fully extended on the floor cotton side up
- Lay patient on the Evac Sledge.
- Ensure patients feet are in the foot pocket.
- Tighten security belts over patient.
- Pull the patient to safety.