

USER INSTRUCTIONS THERAPIST KNEELER

L34277

Therapist Kneeler

⚠ General Safety Warning ⚠

Please make sure that you read and follow these instructions carefully. Keep these instructions in a safe place for future reference.

1. INTENDED USE

Designed to give cushioned comfort to the Therapist whilst kneeling.

A visco elastic memory foam topping ensures ultimate comfort and relief. This lightweight kneeler comes complete with a handy strap for ease of carrying.



2. SPECIFICATION

- Firm Polyurethane foam base, fused with a visco elastic memory foam top.
- Size: 260 x 380 x 40 mm (10 x 15 x 1½").
- Manufactured in UK.
- CE Marked - Fire retardant to source (crib)5.
- 2 way stretch vapour permeable cover.
- Infection Control welded seams with skirted zip.

3. UNPACKING

Remove the packaging and check the product carefully. Contact your supplier immediately if parts are damaged and do not try and use.

4. IN USE

When using the kneeler place the firm foam towards the floor so the knee cushions into the visco topping.

5. CLEANING INSTRUCTIONS.

- Remove foam insert/filling where possible
- Machine wash and tumble dry at 70°C Maximum
- To clean by hand, treat with methylated or white spirit followed by hot water and detergent.

6. PRODUCT DISPOSAL

This product must be thoroughly cleaned/sterilised prior to disposal.

This product is manufactured from foam and fabric and can be disposed of through Local Authority Household Waste and Recycling Centres or through Local Trade Waste Collectors.

7. PRODUCT GUARANTEE

NRS Healthcare guarantee this product for a period of 24 months, from the date of purchase, against faulty materials and workmanship.

This guarantee does not cover abuse, interference or tampering with any part of the products. Issues arising from fair wear & tear and customer misuse are not covered.

This guarantee covers domestic use only. If the product is used in commercial or semi-commercial settings such as care homes, schools, rental or communal accommodation etc, the guarantee is valid for one year only. Consumable parts are also guaranteed for one year only.

Any guarantee issued by NRS Healthcare does not affect your statutory rights.

NRS HEALTHCARE Sherwood House, Cartwright Way, Forest Business Park
Bardon Hill, Coalville, Leicestershire, LE67 1UB, UK.

Tel: Customer Services: +44 (0)345 121 8111

Email: customerservice@nrshealthcare.co.uk

NRS Healthcare is a leading supplier of community equipment and daily living aids.
Please visit one of our websites to view our full range of products:
www.healthcarepro.co.uk or www.completecareshop.co.uk