CT-2 <u>Chair se</u>nsor



CT-2 Easylink chair occupancy detection pressure sensor

Can help reduce the risk of falls when a person rises from their chair (required alarm device)

Position under chair cushion to detect a seated person

Compatible with Easylink-Medpage chair alarm products and transmitters

CT-2 CHAIR SENSOR, corded fall prevention chair occupancy sensor for use with CTM-3 voice alarm transmitter and Easylink-Medpage alarm transmitter. Manufactured in high quality plastics with silver coated internal switch contacts for optimum performance and service life.

Technical Information

Technical Information
CT-2 Chair sensor (380mm x 300mm)
3M connector lead
Connector: 2.5mm mono Jack plug
Silver lined anti-corrosion switch contact
Durable, flexible material
Resistant to diluted disinfectants