



## Boost Adjustable Transfer Aid



**NEW!**  
increased user  
weight capacity of  
200kg

# Boost Adjustable Transfer Aid

Featuring an upgraded leg opening mechanism, increased weight capacity and manoeuvrability enhancements, our NEW Boost Adjustable Transfer Aid truly makes light work of sit to stand transfers!

Specifically designed for those who can raise themselves from a seated position to a standing position, the Boost Adjustable Transfer Aid is ideal for transferring from room to room, enabling users to participate in transfers and maintain patient independence.





## NEW BOOST ADJUSTABLE TRANSFER AID FEATURES!




**NEW! height adjustable leg pad**



**NEW! additional carer push bar**



**NEW! improved opening mechanism**



**NEW! enhanced castor design**

### Boost Adjustable Transfer Aid

**PR31PT007**

- » Lightweight and effortlessly manoeuvrable
- » Tool-free height adjustment of padded lower leg support to suit user comfort
- » Generous seat area for user comfort
- » Upgraded leg-opening mechanism for easy positioning around wider furniture
- » Improved castors for enhanced and secure transportation
- » The Boost can be wheeled directly up to a toilet, with wheels positioned either side of the toilet offering the user more dignity to sit down
- » Narrow width enables the Carer to manoeuvre easily through standard doorways and between furniture
- » Padded hand-rail for user comfort
- » Additional Carer push bar for improved manoeuvring and varying support requirements
- » Easy wipe-clean surfaces
- » Substantial footplate provides stability
- » Robust pedals allow legs to be easily opened and closed smoothly
- » Compatible with our range of optional ladder style support belts for added security

**NEW!**  
increased user  
weight capacity of  
200kg

### PRODUCT SPECIFICATION

Height	<b>104.5cm</b>
Overall Width (legs closed)	<b>66cm</b>
Width (legs open)	<b>91cm</b>
Length	<b>97cm</b>
Seat Height	<b>66cm front, 79.5cm rear</b>
Seat Depth	<b>30.5cm</b>
Seat Width	<b>57cm</b>
Product Weight	<b>33.9kg</b>
User Weight Capacity	<b>200kg (31.5st)</b>

### ACCESSORIES



For additional safety and peace of mind a support belt can be used in combination with the Boost Adjustable Transfer Aid.

The ladder style belts can be locked in place over the belt attachment points for greater stability.