

Sun-Shine Land International Ltd.

8F, No. 69-10, Sec. 2, Chung Cheng East Road, Tamsui Dist., New Taipei City, Taiwan

- T +886-2-2809-5156
- F +886-2-2809-1548
- E info@sunshineland-adl.com

User Instruction for HA-4140 (Economy Week-Day Pill Dispenser)

Please read and retain this instruction, it provides essential information for correct and safe use. This instruction is guideline only.

INTRODUCTION

The economy week-day pill dispenser is ideal for people who take larger numbers of tablets more frequently each day. It has 7 pillboxes each with 4 doses/medication removable compartments per day. Organize all your daily medications once a week. Pill organizer system is great for travel. Helps pill identification and medication adherence. Avoid medication errors.

SPECIFICATION

Overall Size	Weight	Material
226 x 125 x 25 mm	161 g	PP

INSTRUCTIONS FOR USE

- 1. To remove the individual tablet containers, locate the hole underneath and push through with your finger.
- 2. To replace the individual tablet containers, push it back into the cavity.

WARNING

- Keep out of the reach of children.
- Periodically inspect for wear and damage.

MAINTENANCE

The economy week-day pill dispenser can be cleaned by water-washed.

