

P42586



Wheelchair Apron

Tablier de fauteuil roulant

Providing a waterproof covering for the lower body.

Light weight and compact

- Easy to use with an elasticated waist and base
- Attaches easily to the wheelchair with adjustable straps
- Fits over the armrests or around the torso



Fits most manual wheelchairs

Fitting Instructions:

Slip feet in to the base of the Apron and fit around torso.

If required fix to the wheelchair with the adjustable straps.

Care Instructions: -Machine washable at 30°C

-Do not tumble dry - Do not Bleach - Do not iron - Line Dry

Made in China for:

NRS Healthcare,

Sherwood House, Cartwright Way, Forest Business Park, Bardon Hill, Coalville, Leicestershire. LE67 1UB United Kingdom.

www.nrs-uk.co.uk Tel: 0345 121 8111 Fax: 0845 121 8112

