


Note! Read the manual for QuickMove before use.

Adjust the length of the BackSupport to fit the patient. Use of the strap should give the patient support and a sense of security.

Max: 170 kg

Manual no: 00760 En Ver. 1 110506

Assembling QuickMove:




Care of the product:

Read the product label.



Do not use rinsing agent.
To prolong product life, avoid
tumble-drying.



Remove the clasps before laundering and replace them, as illustrated, after laundering.

Contact your local distributor if you have any questions about the product and its use. See www.handicare.com for a complete list of distributors. Always make sure that you have the right version of the manual. The most recent editions of manuals are available for downloading from our website, www.handicare.com.

Handicare AB
Veddestav. 15. Box 640 Tel: +46 (0)8 -557 62 200
SE-175 27 Järfälla Fax: +46 (0)8-557 62 299
SWEDEN

Internet: www.handicare.com
E-mail: info@handicare.se