Good Grips[®] Coated Tablespoon



1. Hold the handle firmly.

The unique Good Grips® utensils have a special twist built into the metal shaft to allow the utensils to be bent to any angle. A great solution for persons with Parkinson's disease, Cerebral Palsy or other neurological impairments.

Not recommended for those with heavy biting reflexes.

Cold-water washing will prolong product life.



2. Bend at the twist, either left or right.



North Coast Medical, Inc.



3. Adapted for custom use.



Printed on Recycled Paper

North Coast Medical, Inc.
San Jose, CA 95125

Good Grips is a registered trademark
of OXO International, handle has a US Patent

GOOD GRIPS®-

A great solution for persons with arthritis, stroke or other neurological impairments.

Good Grips[®] utensils are ideal when bending the wrist, grasping, or other eating motions become difficult. The large handle provides a comfortable, secure grip. The unique twist in the metal shaft allows the utensil to be bent at any angle to accommodate limited wrist motion.

Dishwasher safe



1. Hold the handle firmly.



2. Bend at the twist, either right or left.



-3. Adapted for custom use.

NC 65593 Good Grips® Teaspooth



North Coast

North Coast Medical, Inc. Gilroy, California 95020-7304 U.S.A. www.ncmedical.com

OXO® and Good Grine® and

