

Product Technical Data Sheet 006

Malvern Vinyl Perching Stool Range

Product Range

F19739	Adjustable Height Stool
F19741	Adjustable Height Stool with Armrests
F19893	Adjustable Height Stool with Armrests & Backrest
L33236	Adjustable Height Stool with Armrests & Padded Backrest
F21541	Adjustable Height Stool with Padded Armrests & Padded Backrest
L63344	Extra Low Adj. Height Stool with Armrests & Padded Backrest
L63356	Extra Low Adj. Height Stool with Padded Armrests & Padded Backrest

Product Description & Intended Purpose

These stools are specifically designed to allow users to sit at a comfortable height in the kitchen, bathroom or shower. The padded vinyl seat provides greater comfort and is angled to enable the user to sit and rise with greater ease. The stool is also available with padded vinyl arm rests and back support.

Product Features

- Spayed legs fitted with rubber ferrules provide stability on all floor surfaces.

Dependent on Selected Model

- Padded vinyl seat, arm rests and back support to provide greater comfort.
- Height adjustable mild steel frame.

Technical Specification

- Seat Height: 685mm max x 535mm min
- Width between arms, top: 470mm
- Width between arms, bottom: 435mm
- Backrest height from seat: 270 to base of tube
- Seat: 335 x 290mm
- Backrest: 285 X 140mm
- Armrest Height from seat: 197mm
- Footprint width & depth max: 640W x 460D
- Footprint width & depth min: 560W x 420D
- **Maximum User Weight:** 25st/160kg
- **Product Weight:** 4.5kg

Installation Information *(Note: for detailed information refer to User Instructions)*

- For height adjustment, remove 'e-clips', slide tube to desired position, align holes and re-insert 'e-clips'.

PTDS006

Malvern Vinyl Perching Stool Range

Safety Precautions

- Do not move the stool whilst occupied.
- Do not exceed the maximum user weight of 25st 160kg.
- Ensure the floor surface is level to keep the stool stable and avoid rocking.
- For adjustable models, avoid extending the legs in a manner which causes the stool to slant and ensure legs are fully secured before using the stool.

Cleaning

- To clean the seat, arm and back rests use a mild disinfectant and warm water, wipe seat and dry with a disposable paper towel.
- To clean the steel frame, use a mild disinfectant and warm water.
- Please check the suitability of the cleaning product before use. Do not use abrasive cleaners or cloths as this will permanently damage the surface and invalidate any warranty.
- Product can be autoclaved up to a temperature of 80°C.
- Always allow product to fully dry prior to stacking and storage.

Care & Maintenance

- Check that all fixing clips are inserted correctly and there is no excessive movement.
- Periodically check the frame of the product for signs of wear. Any damage should immediately be reported to the supplier and the product taken out of action.

Accessories / Related Products

None

FAQ's

- How do I know if a perching stool is suitable for me?
The perching stool is designed for use from a semi standing position so you must be able to take some weight through your legs.
- Can the stool be used in the shower?
Yes, the stool has been designed for use in bathrooms, shower rooms and kitchens.

Product Supplied By:

NRS Healthcare Sherwood House, Cartwright Way, Forest Business Park, Bardon Hill, Coalville, Leicestershire. LE67 1UB UK