

Essential Care

Quality Foam Mattress




A cost effective solution for comfort
and pressure redistribution



The **Essential Care** mattress is a pressure redistributing support surface for the community, residential and community healthcare environment. Its castellated surface provides a cost effective solution for comfort and pressure redistribution to meet the needs of individuals at risk or high risk of developing pressure ulcers. The Invacare **Essential Care** is suitable for those considered at 'High Risk' of developing pressure injury.

- ▶ Cut foam castellations
- ▶ Multi-stretch polyurethane fabric cover
- ▶ Zipped around 3 sides for ease of inspection and laundering



Cut out castellations

- ▶ Castellations improve comfort and ensure effective pressure redistribution.

Cover

- ▶ The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability, which helps to reduce shear and friction forces. The cover also has a white inner liner allowing easy inspection for strikethrough or other cover damage.

Design

- ▶ The **Essential Care** mattress has an innovative and modern cover design with the simple wave giving the mattress a fresh and current feel.

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Size range ▼



Max. user weight ▼



Max. cleaning temperature ▼

Nominal density range (kg/m³) ▼

Article no. ▼

Essential Care (150 mm height)	2000 x 880 x 150 mm 1950 x 880 x 150 mm	140 kg	95°	32 - 34	1586659 1626731
Essential Care (120 mm height)	1950 x 880 x 120 mm	110 kg	95°	32 - 34	1626721

Mattresses is manufactured to:
EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com

1. *Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).*